Processing instructions

for the Ligamys suturing forceps

Preservation in motion

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English Page 21–36

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Nederlands Pagina 105–121

Türkçe Sayfa 122–137

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1.	Scope	22
2.	Terminology (EN ISO 17664)	22
2.1	Processing/Reprocessing	22
2.2	Processor	22
2.3	Treatment chemicals	22
2.4	Certification	22
2.5	Decontamination	22
2.6	Disinfection	22
2.7	Manufacturer	22
2.8	Contaminated	22
2.9	Manual cleaning	22
2.10	Cleaning	22
2.11	Washer/Disinfector	22
2.12	Sterile	22
2.13	Sterilisation	22
2.14	Sterilisation medium	22
2.15	Validation	22
3.	Information to be provided by the manufacturer	23
3.1	Instructions for reprocessing	23
3.2	Limitations and restrictions of reprocessing	23
3.2.1	Warnings and precautions	23
3.2.2	Limitations	24
3.3	Preparation at the place of use	24
3.3.1	Processing during and immediately after use	24
3.4	Preparation before cleaning	24
3.5	Cleaning	26
3.5.1	Instructions for manual pre-cleaning	26
3.5.2	Instructions for automated cleaning and disinfection	28
3.6	Control and maintenance	29
3.7	Function testing	29
3.8	Packaging	32
3.9	Sterilisation	32
3.10	Storage	33
3.11	Efficacy of the treatment process	33
3.12	Service life/Number of processing cycles	33
4.	Symbol	34
4.1	Symbol as used by Mathys Ltd Bettlach	34
5.	Customer information	34
6.	Appendix – quick overview	35
6.1	Manual pre-cleaning	35
6.2	Machine cleaning	36
6.3	Sterilisation process with saturated steam	36

Processing instructions for the Ligamys suturing forceps

1. Scope

The present processing instructions in accordance with the requirements of EN ISO 17664 apply to instruments that are to be reused and therefore need to be reprocessed, as well as to medical devices that are sold non-sterile, but to be used in sterile condition. The Ligamys suturing forceps (re-sterilisable surgical instruments) by Mathys Ltd Bettlach belongs to this group.

The process chemistry parameters as well as the equipment of the following processing instructions are recommendations resulting from the findings of the validation process for instrument processing of Mathys Ltd Bettlach.

The processor is responsible for ensuring that the processing actually performed achieves the desired results with the individual equipment, process chemicals and staff in the processing facility. It is not mandatory to use the same chemicals, parameters and technical equipment as in the validation process of Mathys Ltd Bettlach. Equivalent or alternative products may be used with which a successful validated cleaning and sterilisation process can be demonstrated as a result. If the user has a method that is already established and validated and which is proven to lead to the desired result, there is no need for the user to change the method.

2. Terminology (EN ISO 17664)

2.1 Processing / Reprocessing

Any activity, including cleaning, disinfection and sterilisation, that is necessary to prepare a new or used medical device for its intended use.

2.2 Processor

Organisation that is responsible for carrying out the actions necessary to prepare a new or used medical device for its intended use.

2.3 Treatment chemicals

Formulation of compounds intended for use in reprocessing.

Note

These include washing and rinsing agents, surfactants, detergents, disinfectants, enzymatic cleaners and sterilants.

2.4 Certification

Confirmation by testing and furnishing of objective evidence that the specified requirements are met.

2.5 Decontamination

Removal of infectious contaminations. In the case of objects or surfaces, this is done using physical or chemical means, until the former are no longer capable of transmitting infectious particles. Thus, safe handling of surfaces or objects is ensured.

2.6 Disinfection

Method used to reduce the number of viable micro-organisms on a product to a predetermined level that is suitable for the further handling or use of the product.

2.7 Manufacturer

Organisation that is responsible for the design, production, packaging and labelling of a medical device before it is marketed under its own name, regardless of whether these operations are carried out by this organisation itself or by a third party commissioned accordingly.

2.8 Contaminated

Condition after potential or actual contact with micro-organisms.

2.9 Manual cleaning

Cleaning without use of a washer/disinfector.

2.10 Cleaning

Removal of contamination from an object to an extent that is necessary for further processing or for the intended use.

2.11 Washer / Disinfector

Device for cleaning and disinfection of medical devices and other items used in the medical, dental, pharmaceutical and veterinary field.

2.12 Sterile

Free of any viable micro-organisms.

2.13 Sterilisation

Validated method to kill all types of viable micro-organisms on the product.

2.14 Sterilisation medium

Treatment chemical capable of destroying micro-organisms, including viruses, when used in the correct dilution/dosage and applied for the recommended contact time.

2.15 Validation

Documented procedure for collecting, recording and interpreting the results needed to furnish evidence that

a process consistently yields products that meet the specifications.

3. Information to be provided by the manufacturer

3.1 Instructions for reprocessing

These processing instructions for the Ligamys suturing forceps are based on a validated manual/automated processing procedure. A purely manual or purely automated processing method has not been validated by Mathys Ltd Bettlach and does not lead to sufficient cleaning success.

3.2 Limitations and restrictions of reprocessing

3.2.1 Warnings and precautions

Staff who come into contact with contaminated or potentially contaminated surgical instruments must take the generally accepted precautionary measures. When handling instruments with sharp points or edges, special caution is advised.

Personal protective equipment (gown, mask, goggles, visors, gloves, shoes, overshoes, etc.) are required in order to avoid contact with contaminated or potentially contaminated materials, instruments and products.

Particular care must be taken when handling cutting instruments (reamers, drill bits, rasp, chisels), as they pose a risk of injury to the patients on the one hand and to the staff (operating theatre and CSSD personnel) on the other.

It must be clarified in advance that the patients, as well as the staff (operating theatre and CSSD personnel), do not react with allergic reactions due to material intolerance (various steels and plastic materials) upon direct contact with instruments.

For new instruments delivered to the hospital, Mathys Ltd Bettlach recommends triple cleaning before use in order to build up the protective oxide layer.

In the manual cleaning process, low-foam detergents are to be used in order to ensure visibility of the instruments. In manual cleaning with brushes, it is recommended to keep the instruments always below the surface of the cleaning solution. This ensures that no aerosols form and splashes are avoided that might spread contaminants. To prevent accumulation of detergent residues, the detergents must be completely removed by thorough rinsing of the product surfaces.

Do not place heavy instruments on sensitive instruments.

Do not allow contaminated instruments to dry before reprocessing them. This is important because all the steps described below for cleaning and sterilisation are facilitated by preventing blood, body fluids, bones and tissue residues as well as saline or disinfectant from drying on used instruments.

When loading the cleaning trays and baskets, ensure that instruments cannot damage each other, and that they can be easily rinsed all around.

The chloride and iodide ions contained in detergents and disinfectants can cause corrosion. For this reason, contact of the instruments with such agents must be kept as brief as possible. Rinse the instruments thoroughly with tap water and then deionised water (demineralised water) in order to remove all residues. Never leave wet instruments to stand. Rather, dry them immediately. The condensation moisture resulting from sterilisation can be prevented by extension of the drying phase.

Strongly acidic or alkaline detergents or excessive dosages may corrode and destroy the protective passive layer of the instruments or markings. The concentrations and exposure times recommended by the manufacturers must be observed under all circumstances.

For cleaning reusable instruments, Mathys Ltd Bettlach recommends a combined manual/automated cleaning process with a detergent and enzyme solution at an alkaline pH of < 12. It is extremely important to neutralise the alkaline cleaner completely and rinse it thoroughly off the instruments. In automated cleaning, the information provided by the manufacturers of the equipment and cleaning agents must be followed.

Only instruments by Mathys Ltd Bettlach may be used for the placement of implants by Mathys Ltd Bettlach (see the respective surgical technique); instruments by other legal manufacturers may not.

No additional lettering of any kind may be applied to the instruments.

The instruments are packaged separately and delivered in non-sterile condition. The packaging materials must be disposed of in accordance with the local and country-specific regulations.

If these processing instructions are adhered to, neither damage to nor limitation of the service life of the medical devices are to be expected.

3.2.2 Limitations

Patients considered as risk patients with regard to prion diseases such as transmissible spongiform encephalopathy (TSE), Creutzfeldt-Jakob disease (CJD) and its variants (vCJK), must be operated on with single-use instruments whenever possible.

It is extremely important to neutralise alkaline cleaners completely and rinse them thoroughly off the instruments.

The suturing forceps must be carefully inspected after cleaning. In doing so, the functionality of the forceps must be checked. It must be ensured that the cannula tip of the suturing forceps is sharp for use and has not been damaged by improper handling during cleaning. Also use a PDS II 2-0 thread or equivalent thread of thickness 2-0 to ensure the suture-guiding needle is completely patent. If this is no longer the case, contact your local Mathys partner.

Mathys Ltd Bettlach recommends selecting a cleaning agent with enzymatic additives to decompose blood, body fluids and tissues. Please note that some enzyme solutions are intended specifically for decomposing faecal matter or other organic contaminants, and are therefore not suitable for cleaning surgical instruments.

There are no restrictions for the repeated processing of the instruments by Mathys Ltd Bettlach, since it has a low impact on the service life of the instruments. The service life is generally defined by wear and damage from use.

Use of hard water (°dH value > 14) is to be avoided. For the initial rinsing, soft tap water is suitable. Rinse thoroughly with deionised water (demineralised water) to remove any residues. Tap water often contains high concentrations of minerals (e.g. lime or silica), which can be seen on the instrument surface as stains with sharply defined edges.

Tip

Never leave instruments wet, but dry them instantly.

3.3 Preparations at the place of use

3.3.1 Processing during and immediately after use

First, residues of body fluids and tissues must be removed using a specific nylon brush.

Here it should be noted that the water used does not exceed a temperature of 45°C (113°F), since otherwise the contaminations, such as proteins in the blood and tissues, will be fixated on the instrument.

Saline, blood, body fluids, tissue, bone residues or other

organic particles on instruments must be removed before cleaning as soon as possible in order to avoid drying as well as corrosion.

Tip

Immersion of the used instruments into enzyme solutions or into cold deionised water (demineralised water) facilitates cleaning, particularly in the case of instruments featuring complex shapes and areas difficult to access (e.g. cannulated and tubular designs, etc.). The enzyme solutions degrade the protein-containing substances, preventing blood and protein-containing materials from drying on the instruments.

Saline solutions as well as detergents and disinfectants that contain aldehydes, mercury, active chlorine, chloride, bromine, bromide, iodine or iodide are corrosive and must not be used.

Always follow the manufacturer's instructions for preparation and use of the solutions strictly.

Optimal cleaning is ensured if the instruments are properly cleaned within one hour after use in the CSSD to minimise the risk of drying of substances and materials. If it is not possible to clean the instruments within the specified time, Mathys Ltd Bettlach recommends immersing the instruments into an enzyme solution or deionised water (demineralised water) of room temperature, or to wrap them for up to 6 hours into cloths wetted with enzyme solution or deionised water (demineralised water).

For prevention of contamination, the used instruments must be transported in closed or covered containers to the CSSD.

3.4 Preparation before cleaning

Excessive concentrations of the cleaning agents, as well as strongly acidic and alkaline detergents, can damage the protective oxide layer and cause pitting. When using such agents, the concentration and exposure time recommended by the manufacturers must be strictly adhered to. Mathys Ltd Bettlach recommends use of alkaline cleaning agents with a pH < 12.

The concentration, temperature and exposure time for the cleaning agents recommended by the manufacturers must be observed under all circumstances in order to achieve optimal cleaning.

When using dry detergents in powder form, it must be made sure that these are completely dissolved before use to avoid discolouration or corrosion of the instruments.

Strongly contaminated solutions (blood and/or clouding) should be replaced with freshly prepared cleaning solutions.

If water without detergent additives is used for cleaning, Mathys Ltd Bettlach recommends a water temperature of no more than 45°C (113°F), since otherwise blood contained in contaminations would start to coagulate, lead-

ing to strong fixation of proteins on the instrument, which can then be detached only with great effort in the automated cleaning.

The Ligamys suturing forceps must be treated in a two-stage cleaning process, first manually, followed by automated cleaning, in order to achieve the necessary cleaning results.

Table 1: Overview of the reprocessing according to EN ISO 17664

Process			Reusable surgical instruments
Preparation at the place of use	Condition	Wet	• Recommendation Cold deionised water (demineralised water) or enzyme solution (fluid or wetted cloths) • Max. 6 hours
		Dry	• Max. 1 hour
Decontamination	Preparation	• Selection according to the cleaning and disinfection method	
	Cleaning	Manual	–
		Automated	–
		Ultrasound	+
		Combined manual/automated	+
		Alkaline (pH < 12)	+
		Acidic	–
		Neutral to slightly alkaline (pH 7–9.5)	–
	Rinsing		Final rinse with deionised water (demineralised water)
	Disinfection ¹	Chemically max. 60°C (140°F)	–
		Thermally 90°C (194°F)	+
	Drying	T _{max} /time	115°C (239°F)/15 min
Maintenance	Function check		Mandatory
	Care	Sterilit® I JG600	Mandatory
Sterilisation	Moist heat ²		+
	Ethylene oxide		–
	Formaldehyde		–
	Plasma		–




Legend: + validated method, – non-validated method



¹ Non-automatic disinfection process according to «Guideline for Disinfection and Sterilisation in Healthcare Facilities, 2008»

² Preferred sterilisation method according to EN ISO 17664

3.5 Cleaning

3.5.1 Instructions for manual pre-cleaning

No.	Step	Material / Medium	Images
1	Manual removal of all visible contaminations using a plastic brush* below the water surface until no visible residues are present any more.	• Nylon brush • Tap water (cold)	
2	Flush the suture-guiding cannula through the roller housing of the suture forceps with 50 ml of enzyme-containing cleaning solution. Ensure that liquid visibly flows out of the tip of the cannula.	• 50 ml plastic syringe without Luer lock • 1 % deconex® TWIN PH10 and 1 % deconex® TWIN ZYME (v/v) • Deionised water (demineralised water)	
	If the needle is blocked by tissue fragments, the roller housing must be folded up and the cannula flushed by direct application of a syringe. Ensure again that liquid visibly flows out of the tip of the cannula.		
	If upon rinsing of the suture-guiding cannula with enzymatic cleaning solution or deionised water no or only reduced liquid outflow is seen, neither through the roller housing nor with the roller housing folded up, the suturing forceps may not be used and must be replaced.		

No.	Step	Material / Medium	Images
3	The pre-rinsed suturing forceps are treated with ultrasound in the combined cleaning agent/enzyme solution in an ultrasonic bath (30–60 kHz) for 5 minutes. It is to be ensured that the forceps are completely covered with the cleaning agent.	• Ultrasonic bath • 1 % deconex® TWIN PH10 and 1 % deconex® TWIN ZYME (v/v) in osmosis water	
4	Rinse the suturing forceps well under running tap water.	• Tap water (cold)	
5	Rinse the suture-guiding cannula of the suturing forceps with 50 ml of tap water.	• 50 ml plastic syringe without Luer lock • Tap water (cold)	
6	Rinse the suture-guiding cannula of the suturing forceps with 50 ml of deionised water (demineralised water).	• 50 ml plastic syringe without Luer lock • Deionised water (demineralised water)	
7	The suturing forceps must be checked visually for any residues or damage. If there are still any visible residues, these should be removed using a plastic brush under running water, and the manual cleaning steps 3 to 7 should be repeated.	• Visual control	

* Decontaminate brushes after use and sterilise or dispose of them. Do not use steel brushes.

3.5.2 Instructions for automated cleaning and disinfection

The process was carried out at Mathys Ltd Bettlach using a washer and disinfector (WD) by Miele KG (Miele Professional G7836CD) and a combined manual/automated cleaning process with the detergent and enzyme solution deconex® TWIN PH10 and deconex® TWIN ZYME by Borer Chemie AG.

The cleaning basket of the WD must be equipped with a Luer lock attachment and a quiver nozzle.

For correct machine cleaning, the Ligamys suture forceps must be positioned with the front end in the quiver nozzle of the WD to actively clean the suture-guiding cannula. Furthermore, the Ligamys suture forceps must concomitantly be connected to the Luer lock connection of the WD for optimal cleaning of the suture forceps (Fig. 1).


Fig. 1 Mechanical cleaning via Luer lock adapter and introduction of the front end of the Ligamys suture forceps into a quiver nozzle

No.	Step		Equipment / Medium
1	The Ligamys suture forceps is positioned with the front end in a quiver nozzle of the WD and additionally attached, with the roller housing folded up, to the cleaning basket via Luer lock adapter (see Fig. 1). Make sure the Ligamys suture forceps is not damaged during rotation or by the rotor.		Cleaning basket with quiver nozzle and Luer lock attachment
2	Pre-rinse	Duration: 2 minutes	Cold tap water
3	Cleaning process	Dosage / Temperature: 0.5 % deconex® TWIN PH10 at 35°C (95°F) 0.2 % deconex® TWIN ZYME at 40°C (104°F) Duration / Temperature: 10 minutes at 55°C (131°F)	Enzymatic cleaner 0.5 % deconex® TWIN PH10 and 0.2 % deconex® TWIN ZYME, (v/v) in deionised water (demineralised water)
4	Rinse I	Duration: 2 minutes Temperature: Max. 50°C (122°F)	Tap water
5	Rinse II	Duration: 2 minutes Temperature: Max. 40°C (104°F)	Deionised water (demineralised water)
6	Thermal disinfection¹	Duration: 7 minutes Temperature: 90°C (194°F)	Deionised water (demineralised water)
7	Drying	Duration: 15 minutes Temperature: 115°C (239°F)	–
8	The suturing forceps must be checked visually for any residues or damage. If any residues are visible, the entire manual and automated process must be repeated.		Visual control

¹ Thermal disinfection according to DIN EN ISO 15883-1, Annex B and DIN EN ISO 15883-2

3.6 Control and maintenance

After cleaning, the instruments must be completely dry and free of visible and noticeable residues. Critical areas such as handle structures, long and/or thin boreholes or blind holes, joints and complex structures must be treated with special care. To ensure that all contamination has been removed, it is of crucial importance to inspect each instrument carefully and check it for cleanliness as well as water spots (for example, lime or silicate). If any contamination should be discovered adhering to instruments, the complete manual as well as mechanical cleaning and disinfection process must be repeated immediately.

Once the instrument is visually clean, it must undergo maintenance (see arrows in figures 2). To this end, Mathys recommends use of a care product based on paraffin/white oil that is biocompatible, suitable for steam sterilisation, and vapour-permeable, e.g. Sterilit I JG600 (Aesculap AG). Alternative products must be oil-free, free of silicone oil-containing care products, suitable for steam sterilisation and biocompatible (see the «Red Brochure» by the AKI).

For maintenance, the instruments must be cooled to room temperature since otherwise there would be a risk of metal abrasion. The care product must be manually applied specifically, carefully, and drop by drop to hinged or ball bearings of a snap-fit, rotating or joint mechanism and/or sliding surfaces and then distributed evenly by moving the hinges, joints, snap-fit mechanisms, or sliding surfaces. Excess care product must be removed with a lint-free cloth (the manufacturer's instructions must be observed). «Overspraying» the instruments or immersion baths is not recommended by Mathys. Plastic surfaces must not be treated with care products. Observe the expiry date indicated by the manufacturer of the care products.


Fig. 2 Hinge mechanisms to be treated

3.7 Function testing


The suturing forceps must be checked for damage. In particular, the tip of the cannula must be sharp, and the forceps must be able to be closed and opened without any resistance.


Rotation of the wheel on the roller housing must work properly in both directions.

Also use a PDS II 2-0 thread or equivalent thread of thickness 2-0 to ensure the suture-guiding needle is completely patent. If this is no longer the case, please contact your local Mathys partner.

Damage or wear that might impair the function of the instrument must be reported to your local Mathys Partner. The same will decide on repair and/or exchange of the instruments.

Defects and their causes, as well as correct troubleshooting in case of damage, are shown in the table below.

Defect	Cause	Examination and measure
Bent or broken cannula tips 	• Improper handling • Overpacked trays	Can be successfully bent back • Continued use Not possible to bend back or broken • Do not use instrument any longer and discard
Complete opening of the jaws not possible  Correct Incorrect	• Improper handling • Overload during use • Overpacked trays	• Do not use instrument any longer and discard
Complete closing of the jaws not possible  Correct Incorrect	• Improper handling • Overload during use • Overpacked trays	• Do not use instrument any longer and discard
Tube bent 	• Improper handling • Overload during use • Overpacked trays	Can be successfully bent back • Continued use Not possible to bend back or broken • Do not use instrument any longer and discard
Knurled knob blocked 	• Improper handling • Tissue residues in the roller housing	Successful removal of tissue residues with positive functional check • Continued use Functional check negative • Do not use instrument any longer and discard
Luer connection damaged 	• Improper handling • Overpacked trays • Damage due to being dropped	Connection can no longer be attached • Do not use instrument any longer and discard

Defect	Cause	Examination and measure
Labelling no longer legible 	• Worn off due to use • Fading due to the cleaning process, cleaning agents and sterilisation procedure	Article number, lot and CE marking can still be read • Continued use Article number, batch or CE marking can no longer be read • Do not use instrument any longer and discard
Corrosion 	• Unsuitable cleaning agents • Excessive time spent in cleaning agents and disinfectants • Improper processing	• Do not use instrument any longer and discard
Organic residues 	• Unsuitable cleaning agents and disinfectants • Residues allowed to dry for too long prior to processing • Insufficient brushing, rinsing during manual pre-cleaning	• Repeat the processing If repeat processing is not successful • Do not use instrument any longer and discard
Water spots (lime or silicate deposits) 	• Poor water quality • Rinsing with distilled water not performed • Penetration of cleaning agents containing silicate	• Repeat the processing

Instruments must be replaced if:

1. the surfaces look «chalky».
2. they show any signs of damage (e.g. (hairline) cracks, flaking, deformation, blistering).
3. they have excessive shape changes and/or are visibly warped.
4. the lettering, such as item no. or LOT no., is no longer legible. This likewise applies to surgical instruments that do not comprise any plastic materials and are made of steel only.

For replacement, contact your Mathys Partner.

If stains on the medical devices should be recognisable, their cause must first be ascertained. Thus, coloured spots indicate incompatibility with a process chemical or exceedance of an exposure time. White spots are often residues of lime, process chemicals or salts. Corrosion marks should not be underestimated, and affected instruments should be immediately separated from unaffected ones («flash rust» or «rust bloom»).

As damaged instruments can no longer function properly, all reprocessable instruments must be checked for proper functioning after maintenance/care but before sterilisation (see the «Red Brochure» by the AKI).

Markings on the instruments must be legible. This includes scales indicating angles, for determining implant size, of length and/or depth, and of directions such as «left» and «right». If any scales or other markings should not be legible any more, notify your local Mathys Partner promptly for assessment or for replacement of the instruments, respectively.

Please pay particular attention to the following:

1. The instrumentation must be checked for completeness.
2. The instruments in the tray must be checked for correct arrangement.
3. The instruments must be checked for damage (e.g. (hairline) cracks, deformations, changing gaps between metal and plastics, fractures, corrosion or signs of wear) and damaged surfaces. Damage or wear that might impair the function of the instrument must be reported to your local Mathys Partner. The same will decide on repair or exchange of the instruments or entire instrument trays.
4. The functionality of mobile components (e.g. hinge joints, sliding parts, moving parts, etc.) must be checked in order to ensure that the intended movement sequence can be performed and correctly.

5. Long and narrow instruments must be checked for bending.
6. Instruments that consist of several individual components and have to be assembled for function must be checked for correct assembly and functionality after assembly.
7. Drill bits, reamers, rasps and other cutting instruments must be carefully examined for their cutting edges. It must be ensured that the cutting edges are sharp for use, and that no visible or palpable damage is present. This can be done easily using a 10–12 x magnifier.
8. Instruments that are no longer functional must be returned to Mathys for repair or scarping. Before, the instruments must undergo an entire processing cycle in order to eliminate the risk of infection.

3.8 Packaging

The packaging for sterilisation must be suitable for the sterilisation method with moist heat, i.e. permeability of the packaging for the water vapour must be ensured. Moreover, the packaging provides protection during transport and storage.

Therefore, Mathys Ltd Bettlach recommends double packaging of the instrument tray container.

When using sterilisation wrap, this must be free of detergent residues. Mathys Ltd Bettlach discourages use of recyclable fleece.

3.9 Sterilisation

As the processor of the instruments, the user is responsible for performing a procedure validation of all the above steps in order to ensure complete penetration of the sterilising vapour to the instruments.

Furthermore, the user must take protective measures for sharp or potentially dangerous instruments.

Instructions given by the manufacturer of the sterilising equipment must always be followed. If several sets of instruments are sterilised in one sterilisation cycle, the maximum loading of the device in accordance with the manufacturer's instructions must not be exceeded.

For optimum sterilisation and disinfection, the instrument sets must be duly prepared and packaged into the appropriate trays and containers. Only if this is done can the steam spread and penetrate so as to reach all the surfaces. In the case of steam sterilisation it must be ensured that the product is completely dry after sterilisation.

Steam or moist heat is the recommended sterilisation method for instruments by Mathys Ltd Bettlach (EN ISO 17664). Ethylene oxide, gas plasma and dry heat are not recommended as sterilisation methods for the sterilisation of reusable instruments.

For sterilisation of the instruments, in any case the national recommendations/guidelines are to be complied with.

In the following, the minimum sterilisation parameters performed by Mathys Ltd Bettlach with a sterilisation device (Sterimed FAV6767100S) and validated by microbiological examinations to reach an SAL (sterility assurance level) of 10^{-6} are indicated.

Sterilisation process with saturated steam

Type of cycle	Minimum Temperature in °C ⁵ (°F)	Minimum sterilisation time in minutes	Minimum drying time in minutes	Minimum pressure in mbar ^{6, 7}
Pre-vacuum – pulsating vacuum ¹	134 (273)	18	20	≥ 3042
Pre-vacuum – pulsating vacuum (D) ²	134 (273)	5	20	≥ 3042
Pre-vacuum – pulsating vacuum (GB) ^{3, 4}	134 (273)	3	20	≥ 3042

¹ Recommendation of the World Health Organisation (WHO) for steam sterilisation of instruments with potential TSE/CJD contamination

² Hygiene requirements on the reprocessing of medical devices, Federal Institute for Drugs and Medical Devices, 2012

³ Validated sterilisation process with a minimum sterilisation time of 3 minutes at 134°C (273°F) to achieve a Sterility Assurance Level (SAL) of 10^{-6} in accordance with EN ISO 17665-1

⁴ Validation in the original instruments case

⁵ Maximum temperature 137°C according to SN EN 285

⁶ Pressure during the sterilisation phase at 134°C according to ISO 17665-2

⁷ Pressure during the sterilisation phase at 137°C must be ≥ 3318.5 mbar according to ISO 17665-2

3.10 Storage

The packaged, sterile instruments must be stored protected from dust, insects, vermin and direct sunlight in a dry and cool room. This may be accessible only to relevant staff. The equipment for storage and transportation must be such that any disorder, overloading or falling is avoided. Sterile medical devices must never be stored directly on the floor. The instruments must not be stored near chemicals such as active chlorine that release corrosive vapours. The instruments are used in the sequence of the receipt of the goods, and the sterile packaging of the instruments is inspected meticulously for package integrity before opening.

Each user must, adapted to its validated sterilisation process, specify how long the sterile-packaged instruments may be stored before the next use (ISO 58953-8).

If the packaging or a sterile wrap is torn, punctured, visibly damaged or has become wet, the instrument set must be repackaged and

sterilised. In case of signs of open or damaged lid seals or filters on the sterilisation container, too, the instrument set must be re-sterilised and the sterile filter be replaced. For reusable filters, a careful visual inspection must be carried out.

3.11 Efficacy of the treatment process

The processing method recommended in the present processing instructions has been validated. The results meet the requirements with regard to both the limit value and to the guidance value in terms of protein residues in accordance with the guideline by DGKH, DGSV and AKI for mechanical cleaning and thermal disinfection processes for medical devices (D 2596 F).

3.12 Service life / Number of processing cycles

The service life of the Ligamys suturing forceps is generally defined by wear and damage resulting from use. The user is responsible for testing for optimal functioning (e.g. cutting ability), cleanliness and defects (e.g. corrosion) prior to each use.

If the processing is performed according to the present processing instructions, neither damage to nor limitation of the service life of the medical device are to be expected.

4. Symbol

4.1 Symbol as used by Mathys Ltd Bettlach

Caution

5. Customer information




Mathys Ltd Bettlach
Robert Mathys Strasse 5
P. O. Box
2544 Bettlach
Switzerland



Tel. +41 32 644 1 644
Fax +41 32 644 1 161

info@mathysmedical.com
www.mathysmedical.com

6. Appendix – quick overview

6.1 Manual pre-cleaning

1		Remove visible contaminants using a nylon brush below the water surface until no visible residues are present any longer.
2		Flush the suture-guiding cannula with 50 ml of enzyme-containing cleaning solution. Ensure that liquid visibly flows out of the tip of the cannula.
3		If the needle is blocked by tissue fragments, the roller housing must be folded up and the cannula flushed by direct application of a syringe. Ensure again that liquid visibly flows out of the tip of the cannula.
4		If during flushing no liquid should flow out of the front end of suture-guiding cannula, the Ligamys suture forceps may not be used and must be replaced.
5		Insert instruments for at least 5 minutes into an ultrasonic cleaning device (30–60 kHz). Max. temperature 40°C (104°F).
6		Rinse the suturing forceps well under running tap water.

7		Rinse the suture-guiding cannula of the suturing forceps first with 50 ml of tap water. Then rinse with 50 ml of DI water.
8		Check suture forceps visually for any residues or damage. If there are still any visible residues, • remove them using a nylon brush under running tap water and • repeat the entire manual pre-cleaning.

6.2 Machine cleaning

The Ligamys suture forceps is positioned with the front end in a quiver nozzle of the WD and additionally attached, with the roller housing folded up, to the cleaning basket via Luer lock adapter (see Fig. 1). Make sure the Ligamys suture forceps is not damaged during rotation or by the rotor.		• Cleaning basket with quiver nozzle and Luer lock attachment
Pre-rinse	Duration: 2 minutes	• Cold tap water
Cleaning process	Duration: 10 minutes Temperature: At 55°C (131°F)	• Enzymatic cleaner 0.5 % deconex® TWIN PH10 and 0.2 % deconex® TWIN ZYME, (v/v) in deionised water (demineralised water)
Rinse I	Duration: 2 minutes Temperature: Max. 50°C (122°F)	• Tap water
Rinse II	Duration: 2 minutes Temperature: Max. 40°C (104°F)	• Deionised water (demineralised water)
Thermal disinfection	Duration: 7 minutes Temperature: 90°C (194°F)	• Deionised water (demineralised water)
Drying	Duration: 15 minutes Temperature: 115°C (239°F)	–
The suturing forceps must be checked visually for any residues or damage. If any residues are visible, the entire manual and automated process must be repeated.		• Visual control

6.3 Sterilisation process with saturated steam

Type of cycle	Minimum Temperature in °C (°F)	Minimum sterilisation time in minutes	Minimum drying time in minutes	Minimum pressure in mbar
Pre-vacuum – pulsating vacuum ¹	134 (273)	18	20	≥3042
Pre-vacuum – pulsating vacuum (D)	134 (273)	5	20	≥3042
Pre-vacuum – pulsating vacuum (GB) ²	134 (273)	3	20	≥3042

¹ Recommended sterilisation process, ² validated sterilisation process

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