

Surgical technique / Product information

Centris MIS

Preservation in motion

Building on our heritage Moving technology forward Step by step with our clinical partners Towards a goal of preserving mobility

# Preservation in motion

As a Swiss company, Mathys is committed to this guiding principle and pursues a product portfolio with the goal of further developing traditional philosophies with respect to materials or design in order to address existing clinical challenges. This is reflected in our imagery: traditional Swiss activities in conjunction with continuously evolving sporting equipment.

### Table of contents

Inti	roduction	4
1.	Indications and contraindications	5
2.	Preoperative planning	6
3.	Surgical technique	10
4.	Implants	27
<b>5.</b> 5.1	Instruments Additional instruments Centris long	<b>34</b> 38
6.	Measuring templates	40
7.	References	40
8.	Symbols	41

#### Remark

Please make yourself familiar with the handling of the instruments, the productrelated surgical technique and the warnings, the safety notes as well as the recommendations of the instruction leaflet before using an implant manufactured by Mathys Ltd Bettlach. Make use of the Mathys user training and proceed according to the recommended surgical technique.

### Introduction

Nowadays, implantation of an artificial hip joint is considered a routine operation. Hip surgery has three major purposes: relieving pain in the affected hip, restoring the patient's anatomy and function and, improving the range of hip movement. Due to increased life expectancy and broadening of the operative indications, the number of total or partial hip replacements is steadily increasing worldwide.

To prevent complications, a standardized, reproducible and reliable surgical technique is mandatory. The surgical technique provides a stepwise approach for planning and implantation of the Centris Total Hip Arthroplasty.

#### **Centris stem**

The Centris stem in combination with a prosthetic head and acetabular cup component (cemented/uncemented) or native acetabulum represents a system for hip arthroplasty, that is intended to restore the function of the hip joint and/or relieve pain in skeletally mature patients. The Centris stem features a highly-polished stainless steel surface, a rectangular cross section, a CCD angle of 130° and a 12/14 taper.

#### **Charnley-Kerboull Philosophy**

The cemented femoral system is based on the Charnley-Kerboull philosophy. Bone fixation is based on a clinically proven «canal filling stem» concept <sup>1</sup>, the stem being undersized by 0.7 mm compared to the broach. Thus, the stem fills the medullary canal to a large extent and aligns and stabilizes itself during insertion.

#### Design features and advantages of the Charnley-Kerboull Philosophy

- The rectangular cross-section provides rotational stability<sup>1</sup>
- The rounded edges avoid stress concentration in the corners of the cement mantle <sup>1</sup>
- The double-tapered conical shape of the mirror-finished stem converts shear forces into compression forces. Thus, deleterious tensile and bending forces at the stem-cement and bone-cement interfaces can be avoided, ensuring stable long-term fixation of the implant <sup>1</sup>
- The highly polished surface finish with low surface roughness reduces the risk of cracking in the cement mantle <sup>1</sup>

## 1. Indications and contraindications

#### Indications

- Primary or secondary osteoarthritis of the hip
- Necrosis of the femoral head
- Femoral head and femoral neck fractures
- Revision surgery

#### Contraindications

- Presence of factors jeopardizing stable anchoring of the implant:
  - Bone loss and/or bone defects
  - Insufficient bone substance
  - Medullary canal not suitable for the implant
- Local and/or general infection
- Hypersensitivity to any of the materials used
- Severe soft tissue, nerve or vessel insufficiency that jeopardise the function and long-term stability of the implant
- Patients for whom a different type of reconstruction surgery or treatment is likely to be successful

For further information, please refer to the instructions for use or ask your Mathys representative.

## 2. Preoperative planning

Preoperative templating can be performed using conventional radiographs or a digital planning system. The main goal is to plan the appropriate implant size and position to restore the individual biomechanics of the hip joint. Thus, potential problems can be anticipated even before surgery. In most cases, restoration of hip biomechanics can be achieved by reconstructing the original hip rotation centre, the leg length and the femoral and acetabular offset<sup>2</sup>.

Furthermore, the preoperative planning serves as a template in the context of intraoperative balancing by means of fluoroscopic monitoring<sup>3</sup>.

It is recommended to document the preoperative planning in the patient's file.



Fig. 1





Hip templating can best be performed on a pelvic radiograph taken with the patient in standing position. The radiograph needs to be symmetrical, centred on the symphysis of the pubis and with both femora in about 20° of internal rotation. The magnification scale of the radiograph can be controlled with a calibration object or by using a fixed film-to-focus distance and positioning the patient at a fixed distance between film and X-ray source (Fig. 1).

#### Remark

When the affected hip is severely damaged, templating on the unaffected side and mirroring the planning to the affected side should be considered.

The rotation centres of the healthy (A) and affected (A') hip are defined as the centre of a circle that fits the respective femoral head or the acetabular cavity. A first, horizontal line is drawn tangent to both ischial tuberosities, and a second, perpendicular line is plotted through the centre of the symphysis of the pubis.

#### Remark

In case of leg length correction, the adjustment of the leg length can already be considered now, using the ischial tuberosities as a reference.



The acetabular offset can be defined as the distance between Köhler's teardrop (B or B') and a vertical line through the hip rotation centre (A or A') and parallel to the symphysis line (Fig. 2).

#### Planning of the cup

The cup position in relation to the pelvis must take into account the acetabular contours, the hip rotation centre, Köhler's teardrop and the required cup inclination angle (Fig. 3).

Fig. 3



To find an appropriate cup size, various cup templates are positioned at the level of the acetabular cavity aiming to restore the native hip rotation centre while establishing sufficient bone contact, both at the level of the acetabular roof and at the level of Köhler's teardrop (Fig. 4).

Fig. 4



The cup is positioned into the acetabulum. The implant position is established in relation to the anatomical landmarks (acetabular roof, Köhler's teardrop) and the implantation depth is marked down (Fig. 5).

Fig. 5



#### Estimation of the femoral offset

The femoral offset is defined as the smallest distance between the central longitudinal axis of the femur and the hip rotation centre (Fig. 6).

Fig. 6



Fig. 7

#### Planning of the Centris stem

The complete Centris stem system is available in 3 versions Dysplasia (D = 5 sizes), Standard (S = 18 sizes), Revision (R = 4 sizes). Within the Standard version and Revision version two different types of implants are available – a standard type and long stem type. The stem offset is proportional to the stem size (diameter), and various neck length are available for each stem size. The first digit (1 – 5) in the stem identification code indicates the offset (neck length). The middle letter (S, D, R) indicates the stem version, and the last digit (1 – 3) represents the implant diameter.

First, a stem template that restores the original femoral offset is chosen. The stems aligning along the femoral axis and templates with increasing stem offsets (neck length – first digit: 1 to 5) are superimposed on the pelvic radiograph until a match with the original femoral offset is found. Then, within the chosen offset range (first digit), a stem version (second letter: S, D, R) and size (1 to 3) that fills up the femoral canal is selected (Fig. 7).





The correct stem diameter is the diameter that fills the femoral cavity to within 1-2 mm distance from the inner femoral cortex (Fig. 8).

#### Remark

In most cases, an adequate match can be found, i.e. a stem that fills the femoral cavity within 1-2 mm of the cortex and restores both femoral offset and leg length.

Where this is not the case, a compromise needs to be found, or a different hip system should be considered.

Once the final stem has been chosen, the femoral resection level and the stem insertion depth are marked down. To reproduce the stem insertion depth during surgery, the relation of the femoral neck cut to the lesser trochanter, the greater trochanter and the junction between the femoral neck and greater trochanter are registered.

## 3. Surgical technique

Depending on the positioning of the patient and the selection of the approach route, conventional approaches are differentiated from minimally invasive approaches that strive to minimise bone and soft-tissue damage. The Centris stem can be implanted through conventional as well as through minimally invasive approaches. The choice of a specific approach should be based on the patient's anatomy, the personal experience and preferences of the operating surgeon.



Fig. 9



Fig. 10

#### Femoral osteotomy

The femoral neck resection level is related to the distance between the lesser and the greater trochanter and marked according to the preoperative planning (Fig. 9).

#### Remark

The level of the femoral osteotomy be performed proximally to the planned insertion depth of the stem. This allows the femoral neck cut to be finalized with the calcar reamer.

#### Remark

When anatomical conditions prevent head removal after a single neck cut, it is advisable to perform a double osteotomy and remove a fragment of the femoral neck first. Then the femoral head is removed with a femoral head extractor.

Depending on the preference of the surgeon the preparation of the acetabulum and implantation of the cup are to be performed (Fig. 10).

#### Remark

The implantation of the cup is described in a separate surgical technique which can be downloaded from the Mathys Ltd Bettlach website or requested from your local Mathys representative.



Fig. 11



Fig. 12



#### Preparation of the femoral canal

Orthograde implantation is possible only after sufficient lateral opening of the femoral canal. Therefore, the box chisel (Figs. 11-12) must be applied slightly medially of the piriform fossa and introduced in a parallel direction to the dorsolateral femoral cortex with careful hammer strokes.

The opening of the femoral canal with a box chisel should be done carefully, so that there is no fracture of the greater trochanter (Fig. 11).

#### Remark

Pay attention to the desired anteversion of the stem of approximately  $10^{\circ}$ -  $15^{\circ}$  during this step.

The box chisel should be introduced only 1-2 cm proximally into the medullary cavity, otherwise there is a risk of perforation (Fig. 12).

If in doubt, a sharp spoon may be used to explore the inner lateral femoral cortex before the use of the box chisel. Thus, the risk of varus or valgus malposition of the implant is reduced.

Further opening with the reamer facilitates insertion and centring of the subsequent rasps or reamers (Fig. 13).

It must be ensured that the reamer remains its central position aligned to the femoral axis along the inner cortex of the femur as a guide element for preparation of the orthograde reaming.

#### Remark

The cancellous bone must not be removed completely in the process.

Fig. 13



Locking and securing of the smallest rasp in the rasp handle (Fig. 14).

Fig. 14



Fig. 15

Stepwise rasping of the femur.

#### Remark

It is recommended to start with the smallest rasp and then gradually open the femoral canal up to the preoperatively planned size (Fig. 15).

The rasps are introduced along the lateral cortex with moderate hammer strokes into the femoral canal.

#### Remark

The direction of advancement of the rasp needs to be in line with the femoral axis in order to reduce the risk of undersizing or malalignment of the final implant.

Diameter						
	1D1	-	152	-	-	-
	2D1	251	252	253	-	-
Offset	3D1	351	352	353	-	-
	4D1	4S1	452	453	-	-
	5D1	-	-	-	5R1	5R2

Fig. 16

The Centris stem system is available in Standard, Dysplasia and Revision version. Within the Standard version and Revision version two different types of implants are available – a standard type and long stem type (more information about the long stem type can be found on page 17). They are marked as follows:

Standard version	10 sizes (1S2–4S3)
Dysplasia version	5 sizes (1D1–5D1)
Revision version	2 sizes (5R1–5R2)

The first digit (1-5) in the stem identification code indicates the neck length. The middle letter (S, D, R) indicates the stem version, and the last digit (1-3) represents the implant diameter (Fig 16).

During the progressive rasping, select the desired stem version (Standard, Dysplasia, Revision), and increase first the rasp offset up to the planned neck length (first digit). Once the desired offset (neck length) is reached, increase the rasp size (last digit) up to the planned stem diameter.

Do not use rasps of size 3 unless this size has been preoperatively planned as the final implant size. In case of using a Centris long stem, do not use rasps of size 2 unless this size has been preoperatively planned.

#### Example

Planned stem: 4S2 Rasping sequence:  $1S2 \rightarrow 2S1 \rightarrow 3S1 \rightarrow 4S1 \rightarrow 4S2$ 





While progressively widening the medullary canal using increasing rasp sizes, make sure to advance the rasps along the axis of the proximal femur, and control the anteversion of the stem (Fig. 17).

#### Remark

Each rasp should be fully introduced down to the level of the resection plane in order to prevent protrusion of the final implant.



Fig. 18

Fig. 17

Once the largest possible rasp has been introduced down to the femoral resection level or a few millimetres more distally as planned, the connection to the rasp handle is released (Fig. 18).

#### Remark

As soon as cortical contact is perceived rasping must be stopped in order to prevent possible fissures.

If the largest possible rasp is smaller than the stem size that has been templated, early locking of the rasp can be due to:

- 1) Incorrect insertion of the rasp, i.e. varus/valgus or rotational misalignment
- 2) High-density cancellous bone commonly found in young patients
- *3) Inaccurate templating or the use of an incorrect radiographic scale factor*

Insertion of a rasp of a size larger than the one that has been templated can be due to:

- 1) A fracture or fissure of the proximal femur
- 2) Inaccurate templating or the use of an incorrect radiographic magnification factor



Fig. 19



Fig. 20

In each of these cases, intraoperative findings should be compared with the preoperative planning to identify the cause of the mismatch. If needed, appropriate measures to correct the cause of the mismatch should be taken.

The size markings of the rasps match the implant sizes.

#### Remark

Correct fit of the rasp in the femur can additionally be checked under image intensification.

With the largest possible rasp in place, the calcar reamer is positioned over the rasp and the neck resection is finalized (Fig. 19).

Thus, the small medial collar of the Centris stem will fit the medial neck cut, allowing adequate control of the insertion depth of the final implant.

The planned and matching trial cone is placed on the rasp (Fig 20).

#### Remark

It is important to ensure the use of the correct trial cone, which corresponds to the rasp family as shown in the table below:

Stem Family	Trial Neck 1 51.34.0917	Trial Neck 2 51.34.0918	Trial Neck 3 51.34.0919	Trial Neck 4 51.34.0920	Trial Neck 5 51.34.0921
1XY	<b>~</b>	-	-	-	-
2XY	-	<b>~</b>	-	-	-
3XY	-	-	<b>~</b>	-	-
4XY	-	-	-	<b>~</b>	-
5XY	-	-	-	-	<b>~</b>



Fig. 21



Fig. 22



The selected trial head, with a diameter corresponding to the inner diameter of the cup, is positioned on the taper (Fig. 21).

#### Remark

*Trial heads for trial reduction are available in diameters* 28 mm, 32 mm and 36 mm, each with neck lengths S, M, L, XL and XXL.

An overview of the neck lengths of the trial heads can be found in the chapter «Instruments».

Before trial reduction, it is recommended to compare the position of the rotation centre of the trial head and the insertion depth of the rasp with the preoperative templating measurements.

#### Remark

Match the final head size with the inner diameter of the cup.

Trial reduction with the final rasp (Fig. 22).

After trial reduction, take the hip through a full range of motion. Watch out for soft tissue and neck-cup impingement and evaluate the tendency of the implant to dislocate during internal and external rotation in flexion and extension. Also make sure that soft tissue tension is appropriate (Fig. 23).

#### Remark

At this time, it is still possible to modify the stem size and offset, the neck length of the trial head and to some extent the stem anteversion if needed.

#### Remark

Correct fit of the rasp in the femur can additionally be checked under image intensification.

Fig. 23



Fig. 24





Fig. 25

## Additional information for using the Centris long stem

#### Remark

The procedure for implanting a Centris long stem is to start with the rasping process before further distal preparation of the femoral canal with the flexible reamer.

Stepwise rasping of the femur (Fig. 24).

#### Remark

The Centris long stems are available only in the following versions:

	180 mm	230 mm
351	<b>~</b>	<b>~</b>
352	<b>~</b>	<b>~</b>
451	¥	¥
452	<b>~</b>	<b>~</b>

1	175 mm	225 mm
5R1	✓	<b>~</b>

#### Remark

It is recommended to start with the smallest rasp and then gradually open the femoral canal up to the preoperatively planned size.

The rasps are introduced along the lateral cortex into the femoral canal with moderate hammer strokes (Fig 25).

#### Remark

The direction of advancement of the rasp needs to be in line with the femur axis, to reduce the risk of undersizing or malalignment of the final implant.



Fig. 26



Fig. 27

During the progressive rasping, select the desired stem version and increase first the rasp offset up to the planned neck length (first digit). Once the desired offset (neck length) is reached, the rasp size (last digit) is increased up to the planned stem diameter.



Do not use rasps of size 2 unless this size has been preoperatively planned as the final implant size.

#### Example

Planned stem: 4S2 Rasping sequence:  $1S2 \rightarrow 2S1 \rightarrow 3S1 \rightarrow 4S1 \rightarrow 4S2$ 

While progressively widening the medullary canal using increasing rasp sizes, make sure to drive the rasps along the axis of the proximal femur and control the anteversion of the stem (Fig. 26).

#### Remark

Each rasp should be fully introduced down to the level of the resection plane in order to prevent protrusion of the final implant.

With the largest possible rasp in place, the calcar reamer is positioned over the rasp, and the neck resection is finalized (Fig. 27).

Thus, the small medial collar of the Centris long stem will fit the medial neck cut, allowing adequate control of the insertion depth of the final implant.



Fig. 28





Fig. 30

After rasping, use the cylindrical flexible reamers to remove distal cancellous bone only. Stop reaming upon contact with the cortical bone (Figs. 28–30).

Stem size and length		Recommended final flexible reamer size
3S1	180 mm	10
351	230 mm	10
3S2	180 mm	10
3S2	230 mm	10
4S1	180 mm	11
4S1	230 mm	11
4S2	180 mm	11
4S2	230 mm	11

5R1	175 mm	11/12	
5R1	225 mm	11/12	



#### Remark

Mark the length of the trial prosthesis on the flexible reamer in order to check the depth of the implant (Fig 31).





Verification of the correct fit by introduction of the Centris long trial stem (Figs. 32-33).

#### Remark

The size markings of the trial prosthesis match the implant sizes.

#### Remark

Correct fit of the trial prosthesis in the femur can additionally be checked under image intensification.

Fig. 32







Fig. 34



Fig. 35



The selected trial head, with a diameter corresponding to the inner diameter of the cup, is positioned on the taper of the trial stem (Figs. 34-35).

#### Remark

For the final implant, as well as for the trial prosthesis, the standard trial heads with the following item numbers must be used:

Item no.	Description
51.34.1064	Trial head 28 S
51.34.1065	Trial head 28 M
51.34.1066	Trial head 28 L
51.34.1067	Trial head 28 XL
51.34.1068	Trial head 28 XXL
51.34.1069	Trial head 32 S
51.34.1070	Trial head 32 M
51.34.1071	Trial head 32 L
51.34.1072	Trial head 32 XL
51.34.1073	Trial head 32 XXL
51.34.1074	Trial head 36 S
51.34.1075	Trial head 36 M
51.34.1076	Trial head 36 L
51.34.1077	Trial head 36 XL
51.34.1078	Trial head 36 XXL

An overview of the neck lengths of the trial heads can be found in chapter «Instruments».

Before trial reduction, it is recommended to compare the position of the rotation centre of the trial head and the insertion depth of the trial stem with the preoperative templating measurements.

#### Remark

Match the final head size with the inner diameter of the cup.





Fig. 37

#### Trial reduction of trial stem

After trial reduction, take the hip through a full range of motion. Watch out for soft-tissue and neck-cup impingement and evaluate the tendency of the implant to dislocate during internal and external rotation in flexion and extension. Also make sure that soft tissue tension is appropriate (Fig. 36-37).

#### Remark

At this time, it is still possible to modify the stem size and offset, the neck length of the trial head, and to some extent the stem anteversion if needed.

#### Remark

Correct fit of the trial prothesis in the femur can additionally be checked under image intensification.

#### Selection and insertion of the cement restrictor

The inner diameter of the medullary canal can be evaluated preoperatively on radiographs of the proximal femur or using a measuring cone positioner to determine the position and size of the medullary plug.

#### Remark

The measurement is performed on the medial line indicating the resection plane. The medullary plug should be placed 1 cm distally to the tip of the prosthesis.

The medullary plug made of autologous cancellous bone, polyethylene or resorbable synthetic material is used according to the height of the test implantation.

#### Remark

The instruments to be used to determine the size of the medullary plug are not included in the standard instrument set and must therefore be ordered separately.

For more information about the Mathys medullary plug, please contact your local Mathys representative.

#### **Final canal preparation**

Prior to cementing, the medullary canal must be cleaned from loose bone debris and fat that could compromise proper interlocking of the cement with the cancellous bone of the proximal femur. This can be achieved with a curette or brush and extensive jet lavage. It is important to preserve the layer of cancellous bone that is well-fixed to the inner cortex to allow proper cement interlocking during cement pressurisation.

Subsequently the prosthetic bed is carefully aspirated and dried. Concomitantly, the bone cement is mixed.



Fig. 38

#### Retrograde cementing technique

The long nozzle of the cement gun, filled with cement, is introduced against the cement restrictor plug. As cement is injected in a retrograde fashion into the medullary canal, the nozzle is removed progressively until the canal is filled up to the neck resection level (Fig. 38).

#### Antegrade cementing technique

First, a venting tube is introduced against the cement restrictor plug. Then, the short nozzle of the cement gun, filled with cement, is introduced in the proximal femur, and cement is injected into the medullary canal up to the neck resection level. The venting tube prevents air, blood and fat from being trapped by the cement restrictor plug and must be removed before cement pressurization.

#### **Cement pressurization**

To improve cement interlocking, a proximal seal is used to obliterate the proximal femoral canal and additional cement is injected under pressure with the cement gun.

#### Remark

The cementing technique requires special precautions (preparation of the bone cavity, cementing technique, collaboration with the anaesthetist, etc.) which are described in the specific instructions for use of the cement.



Fig. 39



Fig. 40





#### Implantation of the Centris stem

The selected Centris stem, together with the cone protector, is inserted without a hammer into the medullary canal that has been filled with cement.

To avoid misalignment and trapping of air between the stem and the cement, the stem must be introduced slowly and in line with the proximal axis of the femoral canal until the small medial collar reaches the neck resection level (Figs. 39–40).

#### Remark

Alternatively, the stem can be introduced by hand and maintained with the stem introducer.

Remove all excess bone cement and hold the stem steadily in position with the stem introducer or the stem introducer until the bone cement is fully polymerized (Fig. 41).

Remove the cone protection.



Fig. 42







Fig. 44

After curing of the cement, the trial heads of different lengths can be used to perform another trial reduction with the implant in place for testing the range of motion and the ligament tension (Figs. 42-43).

#### Remark

At this time, only the neck length of the prosthetic head can still be modified if needed.

#### Remark

An overview of the neck lengths of the final heads and trial heads can be found in chapter «Implants» and «Instruments».

#### Remark

The head diameter must always match the inner diameter of the cup.

To avoid complications at the stem / head interface, the stem cone needs to be dry and free of any foreign matter (e.g. tissue parts, bone or cement particles) before assembling the final head (Figs. 44–45).



Fig. 45

Reduction of the joint (Figs. 46-47).

#### Remark

Correct fit of the implants can additionally be checked under image intensification.

The joint space needs to be free of any cement or bone particles present.

Depending on the approach, the muscle insertions are reattached, and the wound is closed layer by layer.



Fig. 46

#### Removal of the Centris Stem

In case of revision, the Centris stem can be removed with universal stem extraction instruments and general instruments to remove bone cement.

For further information about stem revision and extraction instruments, please contact your local Mathys representative.



In case of intraoperative removal of the final stem, re-implantation of the same stem is not allowed – a new stem must be used.



Fig. 47

### 4. Implants





ltem no.	Size	Neck length (a)	Stem length (b)
56.11.0055	1D1	28 mm	105 mm
56.11.0056	2D1	32 mm	112 mm
56.11.0057	3D1	35 mm	120 mm
56.11.0058	4D1	39 mm	129mm
56.11.0059	5D1	47 mm	128mm

Material: FeCrNiMnMoNbN Cone: 12/14mm CCD-angle: 130°

#### Centris standard stem, cemented

ltem no.	Size	Neck length (a)	Stem length (b)
56.11.0060	152	28 mm	104 mm
56.11.0061	251	32 mm	112 mm
56.11.0062	252	32 mm	113 mm
56.11.0063	253	32 mm	113 mm
56.11.0064	351	35 mm	121 mm
56.11.0065	352	35 mm	122 mm
56.11.0066	353	35 mm	123 mm
56.11.0067	4S1	39 mm	128 mm
56.11.0068	452	39 mm	129mm
56.11.0069	453	39 mm	129mm

Material: FeCrNiMnMoNbN Cone: 12/14mm CCD-angle: 130°

#### Centris revision stem, cemented

ltem no.	Size	Neck length (a)	Stem length (b)
56.11.0070	5R1	47 mm	131 mm
56.11.0071	5R2	47 mm	133 mm

Material: FeCrNiMnMoNbN Cone: 12/14mm CCD-angle: 130°



#### Centris standard long stem, cemented

ltem no.	Size	Neck length (a)	Stem length (b)
56.11.0072	3S1	35 mm	180 mm
56.11.0074	3S2	35 mm	180 mm
56.11.0076	4S1	39 mm	180 mm
56.11.0078	4S2	39 mm	180 mm
56.11.0073	3S1	35 mm	230 mm
56.11.0075	352	35 mm	230 mm
56.11.0077	4S1	39 mm	230 mm
56.11.0079	4S2	39 mm	230 mm

Material: FeCrNiMnMoNbN Cone: 12/14 mm CCD-angle: 130°

#### Centris revision long stem, cemented

Item no.	Size	Neck length (a)	Stem length (b)
56.11.0080	5R1	47 mm	175 mm
56.11.0081	5R1	47 mm	225 mm

Material: FeCrNiMnMoNbN Cone: 12/14mm CCD-angle: 130°

#### Femoral Heads

#### Femoral Head, Stainless Steel



ltem no.	OD	Neck leng	Jth
54.11.1031	22.2 mm	S	- 3 mm
54.11.1032	22.2 mm	Μ	0 mm
54.11.1033	22.2 mm	L	+ 3 mm
2.30.410	28 mm	S	-4mm
2.30.411	28 mm	Μ	0 mm
2.30.412	28 mm	L	+4mm
2.30.413	28 mm	XL	+8mm
2.30.414	28 mm	XXL	+ 12 mm
2.30.400	32 mm	S	-4mm
2.30.401	32 mm	Μ	0 mm
2.30.402	32 mm	L	+4mm
2.30.403	32 mm	XL	+8mm
2.30.404	32 mm	XXL	+ 12 mm
Marchard La Da Cuntin Asin Asin Ibini			

Material: FeCrNiMnMoNbN Cone: 12/14mm

Femoral Head, CoCrMo



ltem no.	OD	Neck leng	yth
52.34.0125	22.2 mm	S	- 3 mm
52.34.0126	22.2 mm	Μ	0 mm
52.34.0127	22.2 mm	L	+ 3 mm
2.30.010	28mm	S	-4mm
2.30.011	28mm	Μ	0 mm
2.30.012	28mm	L	+4mm
2.30.013	28mm	XL	+8mm
2.30.014	28mm	XXL	+ 12 mm
2.30.020	32 mm	S	-4mm
2.30.021	32 mm	Μ	0 mm
2.30.022	32 mm	L	+4mm
2.30.023	32 mm	XL	+8mm
2.30.024	32 mm	XXL	+ 12 mm
52.34.0686	36 mm	S	-4mm
52.34.0687	36 mm	Μ	0 mm
52.34.0688	36 mm	L	+4mm
52.34.0689	36 mm	XL	+8mm
52.34.0690	36 mm	XXL	+ 12 mm

Material: CoCrMo Cone: 12/14mm

#### **Femoral Heads**

Femoral Head, ceramys



Item no.	OD	Neck leng	th
54.47.0010	28mm	S	-3.5 mm
54.47.0011	28 mm	Μ	0 mm
54.47.0012	28 mm	L	+3.5mm
54.47.0110	32 mm	S	-4mm
54.47.0111	32 mm	Μ	0 mm
54.47.0112	32 mm	L	+4mm
54.47.0113	32 mm	XL	+8mm
54.47.0210	36 mm	S	-4mm
54.47.0211	36 mm	Μ	0 mm
54.47.0212	36 mm	L	+4mm
54.47.0213	36 mm	XL	+8mm

 $\begin{array}{l} \textbf{Material:} \ ZrO_2\text{-}Al_2O_3\\ \textbf{Cone:} \ 12/14\,mm \end{array}$ 

For ceramic-ceramic pairings, use only ceramic heads with ceramic inlays by Mathys.



#### Femoral Head, symarec

Item no.	OD	Neck len	gth
54.48.0010	28 mm	S	-3.5mm
54.48.0011	28 mm	Μ	0 mm
54.48.0012	28 mm	L	+3.5mm
54.48.0110	32 mm	S	-4mm
54.48.0111	32 mm	Μ	0 mm
54.48.0112	32 mm	L	+4mm
54.48.0113	32 mm	XL	+8mm
54.48.0210	36 mm	S	-4mm
54.48.0211	36 mm	Μ	0 mm
54.48.0212	36 mm	L	+4mm
54.48.0213	36 mm	XL	+8mm

**Material:** Al<sub>2</sub>O<sub>3</sub>-ZrO<sub>2</sub> **Cone:** 12/14 mm

For ceramic-ceramic pairings, use only ceramic heads with ceramic inlays by Mathys.

#### **Revision Heads**

**Revision Head, ceramys** 



ltem no.	OD	Neck len	gth
54.47.2010	28 mm	S	-3.5 mm
54.47.2020	28 mm	Μ	0 mm
54.47.2030	28 mm	L	+3.5mm
54.47.2040	28 mm	XL	+7mm
54.47.2110	32 mm	S	-3.5 mm
54.47.2120	32 mm	Μ	0 mm
54.47.2130	32 mm	L	+3.5 mm
54.47.2140	32 mm	XL	+7mm
54.47.2210	36 mm	S	-3.5 mm
54.47.2220	36 mm	Μ	0 mm
54.47.2230	36 mm	L	+3.5 mm
54.47.2240	36 mm	XL	+7mm

**Material:** ZrO<sub>2</sub>-Al<sub>2</sub>O<sub>3</sub>, TiAl6V4 **Cone:** 12/14 mm

ceramys Revision Heads can be used with all Mathys stem systems with a «12/14 cone».

The ceramys Revision Heads can be only combined with Mathys polyethylene or ceramic cups or inlays.

#### Bipolar Head, CoCrMo and Stainless Steel



Material CoCrMo: CoCrMo; UHMWPE Material stainless steel: FeCrNiMnMoNbN; UHMWPE

Detailed information on the implantation of bipolar heads is provided in a separate surgical technique. Please contact your local Mathys agency for this.





#### Hemiprosthesis Head, Stainless Steel Sizes 38–44 mm

Item no. / S -4mm	ltem no. / M 0mm	OD
2.30.420	67092	38 mm
2.30.421	67093	40 mm
2.30.422	67094	42 mm
2.30.423	67095	44 mm

Material: FeCrNiMnMoNbN Cone: 12/14mm



#### Hemiprosthesis Head, Stainless Steel Sizes 46–58mm

Item no. / S -4mm	ltem no. / M 0mm	OD
2.30.424	67096	46 mm
2.30.425	67097	48 mm
2.30.426	67098	50 mm
2.30.427	67099	52 mm
2.30.428	67100	54 mm
2.30.429	67101	56 mm
2 30 430	67102	58 mm

Material: FeCrNiMnMoNbN Cone: 12/14 mm

## 5. Instruments

#### Centris MIS Instrumentation: 51.34.0927A



Item no.	Description
51.34.0924	Centris MIS Tray
51.34.0925	Centris MIS Insert
51.34.0926	Centris MIS Lid





ltem no.	Description
51.34.0900	Centris Rasp MIS 1S2
51.34.0901	Centris Rasp MIS 2S1
51.34.0902	Centris Rasp MIS 2S2
51.34.0903	Centris Rasp MIS 2S3
51.34.0904	Centris Rasp MIS 3S1
51.34.0905	Centris Rasp MIS 3S2
51.34.0906	Centris Rasp MIS 3S3
51.34.0907	Centris Rasp MIS 4S1
51.34.0908	Centris Rasp MIS 4S2
51.34.0909	Centris Rasp MIS 4S3
51.34.0910	Centris Rasp MIS 5R1
51.34.0911	Centris Rasp MIS 5R2
51.34.0912	Centris Rasp MIS 1D1
51.34.0913	Centris Rasp MIS 2D1
51.34.0914	Centris Rasp MIS 3D1
51.34.0915	Centris Rasp MIS 4D1
51.34.0916	Centris Rasp MIS 5D1



ltem no.

51.34.0075

51.34.0076

51.34.0034



51.34.0189	twinSys double offset adaptor right
51.34.0190	twinSys double offset adaptor left
ltem no.	Description
51.34.0033	twinSys calcar reamer 30 mm

twinSys calcar reamer 40 mm

twinSys rasp handle MIS II offset twinSys rasp handle MIS II straight

Description



ltem no.	Description	Neck length
51.34.1064	Trial head 28 S	-4mm
51.34.1065	Trial head 28 M	0 mm
51.34.1066	Trial head 28 L	+4mm
51.34.1067	Trial head 28 XL	+8mm
51.34.1068	Trial head 28 XXL	+ 12 mm
51.34.1069	Trial head 32 S	-4mm
51.34.1070	Trial head 32 M	0 mm
51.34.1071	Trial head 32 L	+4mm
51.34.1072	Trial head 32 XL	+8mm
51.34.1073	Trial head 32 XXL	+ 12 mm
51.34.1074	Trial head 36 S	-4mm
51.34.1075	Trial head 36 M	0 mm
51.34.1076	Trial head 36 L	+ 4 mm
51.34.1077	Trial head 36 XL	+8mm
51.34.1078	Trial head 36 XXL	+ 12 mm

#### **General Instruments**

Item no.	Description
3.30.130	Ruler length 20
ltem no.	Description
3.30.552	Crossbar long
Item no.	Description
3.30.537	Impactor top 36
3.30.538	Impactor top 28
3.30.539	Impactor top 32
ltem no.	Description
51.34.0134	Box chisel silicone
Item no.	Description
51.34.0135	Head impactor silicone
ltem no.	Description
51.34.0263	Impactor/extractor silicone
ltem no.	Description
51.34.0469	Opening reamer for straight stems
Item no.	Description
3.30.536	Top for head impactor
Item no.	Description
56.02.2016	Reamer, narrow
Item no.	Description
58.02.4030	Box chisel MIS



----

Centris MIS – **37** 

### 5. Instruments

### 5.1 Additional instruments Centris long

#### Centris Instrumentation LS Addition 56.01.0013A\*

Tray no.	Item no.
7	56.03.6107
8	56.03.6108
9	51.34.0074

**Note:** The mentioned Item numbers refer to the empty trays. The pictures are intended as support and illustrate the different set options.

\* Optional – needed for long stems



Centris trial pro	sthesis	Tray no.
175 mm	5R1	7
180 mm	351	7
180 mm	352	7
180 mm	4S1	7
180 mm	452	7
225 mm	5R1	8
230 mm	351	8
230 mm	352	8
230 mm	4S1	8
230 mm	452	8
	Centris trial pro 175 mm 180 mm 180 mm 180 mm 225 mm 230 mm 230 mm 230 mm	Centris trial prosthesis   175 mm 5R1   180 mm 3S1   180 mm 4S1   180 mm 4S1   180 mm 4S2   225 mm 5R1   230 mm 3S1   230 mm 4S1   230 mm 4S1   230 mm 4S1   230 mm 4S1

Set no. 56.01.0013A

ltem no.	Reamer flexible 2 <sup>nd</sup> gen.	Tray no.
51.34.0063	8.0 mm	9
51.34.0064	8.5 mm	9
51.34.0065	9.0 mm	9
51.34.0066	9.5 mm	9
51.34.0067	10.0 mm	9
51.34.0068	10.5 mm	9
51.34.0069	11.0 mm	9
51.34.0070	11.5 mm	9
51.34.0071	12.0 mm	9
51.34.0072	12.5 mm	9
51.34.0073	13.0 mm	9
Set no. 56 01 0013A		

Set no. 56.01.0013A

ltem no.	Description
56.02.6183	Guide for medullary reamer flex.
Set no. 56.01.0013A	



ltem no.	Description
58.02.4008	Handle with quick coupling
Set no. 56.01.0013A	

## 6. Measuring templates



Item no.	Description
330.010.007	Centris stem dysplasia template
330.010.008	Centris dysplasia 28 mm template
330.010.010	Centris standard 28 mm template
330.010.011	Centris long stem 28 mm template

### 7. References

- <sup>1</sup> Scheerlinck et al (2006) The design features of cemented femoral hip implant; J Bone Joint Surg [Br] 2006;88-B:1409-18
- <sup>2</sup> Scheerlinck Th. (2010) Primary hip arthroplasty templating on standard radiographs. A stepwise approach; Acta Orthop. Belg., 2010, 76, 432-442
- <sup>3</sup> Loweg L., Kutzner K.P., Trost M., Hechtner M., et al. The learning curve in short-stem THA: influence of the surgeon's experience on intraoperative adjustments due to intraoperative radiography. European Journal of Orthopaedic Surgery & Traumatology, 2017

## 8. Symbols



Centris MIS – **41** 

### Notes




Australia	Mathys Orthopaedics Pty Ltd Lane Cove West, NSW 2066 Tel: +61 2 9417 9200 info.au@mathysmedical.com	Italy	Mathys Ortopedia S.r.l. 20141 Milan Tel: +39 02 5354 2305 info.it@mathysmedical.com
Austria	Mathys Orthopädie GmbH 2351 Wiener Neudorf Tel: +43 2236 860 999 info.at@mathysmedical.com	Japan	Mathys KK Tokyo 108-0075 Tel: +81 3 3474 6900 info.jp@mathysmedical.com
Belgium	Mathys Orthopaedics Belux N.VS.A. 3001 Leuven Tel: +32 16 38 81 20 info.be@mathysmedical.com	New Zealand	Mathys Ltd. Auckland Tel: +64 9 478 39 00 info.nz@mathysmedical.com
France	Mathys Orthopédie S.A.S 63360 Gerzat Tel: +33 4 73 23 95 95 info.fr@mathysmedical.com	Netherlands	Mathys Orthopaedics B.V. 3001 Leuven Tel: +31 88 1300 500 info.nl@mathysmedical.com
Germany	Mathys Orthopädie GmbH «Centre of Excellence Sales» Bochum 44809 Bochum Tel: +49 234 588 59 0 sales.de@mathysmedical.com	P. R. China	Mathys (Shanghai) Medical Device Trading Co., Ltd Shanghai, 200041 Tel: +86 21 6170 2655 info.cn@mathysmedical.com
	«Centre of Excellence Ceramics» Mörsdorf 07646 Mörsdorf/Thür. Tel: +49 364 284 94 0 info.de@mathysmedical.com	Switzerland	Mathys (Schweiz) GmbH 2544 Bettlach Tel: +41 32 644 1 458 info@mathysmedical.com
	«Centre of Excellence Production» Hermsdorf 07629 Hermsdorf Tel: +49 364 284 94 110 info.de@mathysmedical.com	United Kingdom	Mathys Orthopaedics Ltd Alton, Hampshire GU34 2QL Tel: +44 8450 580 938 info.uk@mathysmedical.com

Local Marketing Partners in over 30 countries worldwide ...

Mathys Ltd Bettlach • Robert Mathys Strasse 5 • P.O. Box • 2544 Bettlach • Switzerland

**€€** 0123