Manufacturer

GOMINA AG

Distributed by



For healthcare professional use only. The illustrated image does not represent a con between the use of the medical device described, nor its performance.

Handling brochure

optimys chisel system

Preservation in motion

Building on our heritage Moving technology forward Step by step with our clinical partners Towards a goal of preserving mobility

Preservation in motion

As a Swiss company, Mathys is committed to this guiding principle and pursues a product portfolio with the goal of further developing traditional philosophies with respect to materials or design in order to address existing clinical challenges. This is reflected in our imagery: traditional Swiss activities in conjunction with continuously evolving sporting equipment.

Table of contents

5.	Symbols	24
4. 4.1 4.2 4.3	Instruments Basic instr. set for optimys stem removal 51.34.1173A – Configuration Detail Instrument List (re-usable instruments) Detail Instrument List (single use instruments)	20 21 22 22
3. 3.1	Surgical technique Surgical steps in detail	6 8
2.	Medical Advisor Team	5
1.	Introduction	4

Remark

Please make yourself familiar with the handling of the instruments, the product-related handling brochure and the warnings, the safety notes as well as the recommendations of the instruction leaflet.

1. Introduction

The introduction of bone preserving short stems has created an associated demand for bone preserving revision instruments to remove these stems. To this end our medical advisors, with the close cooperation of the company Gomina, have developed a chisel system specifically designed to extract optimys stems whilst preserving as much bone as possible.



2. Medical Advisor Team



Prim. Priv.-Doz. Dr. Josef Hochreiter (Linz, Austria)



Prof. Dr. Karl Stoffel (Basel, Switzerland)



PD. Dr. Karl Philipp Kutzner (Wiesbaden, Germany)

3. Surgical technique

The optimys stem can be extracted using various surgical approaches and patient positioning. The decision in favour of a specific technique should be based on the patient's anatomy, the reason for revision and the personal experience and preferences of the operating surgeon.



Step 1 Attachment of the Chuck to the cone of the stem



Step 2

Insertion of the Prestarter Chisels through the guide slots of the Chuck (lateral/medial)



Step 3

Insertion of the Starter Chisels through the guide slots of the Chuck (lateral/ventral/dorsal)



Step 4 Removal of the Chuck



Step 5 Application of the Final Chisel

Step 6 Extraction of the stem

3.1 Surgical steps in detail



Layer-by-layer preparation of the surgical field and access to the joint.

Dislocate the joint which has the implanted optimys stem.

Clean the cone and the neck of the stem and expose the lateral portion of the stem shoulder (Fig. 1).



Laterally, the shoulder of the stem must be completely free from bone and soft tissue until the Chuck rests flush on the shoulder and the nose of the Chuck rests in the impaction hole of the optimys stem.

Fig. 1





Place the Chuck that matches the implanted stem (Fig. 2).

Remark

An individual chuck size is required for each stem size and offset version.



Remark

In its final position, the Chuck is flush with the end face of the stem cone (Fig. 3).

Fig. 3



Tighten the fixation screw on the Chuck with the screwdriver (58.02.4005) until the Chuck is firmly fixed to the implant (Fig. 4).

Remark

A firm connection between the Chuck and the implant is necessary to subsequently allow the guided insertion of the curved Prestarter Chisels and Starter Chisels along the curved stem.

Fig. 4



Always start with the lateral Prestarter Chisel.

Insert the Prestarter Chisel into the Handle (Fig. 5).

Always ensure that the chisels are placed in the Handle with the numbers facing upwards (the numbers must be legible when the clamping lever is closed).

Fig. 5



The Handle provides three different clamping positions. The Prestarter Chisels must always be inserted into clamping position 1 at the beginning of the preparation of the respective contour, to ensure sufficient stability of the chisels. The number <1> must be legible through the round borehole of the chisel (Fig. 6).

The chisels must always be inserted flat into the undulating contour of the Handle with the numbers facing upwards.

When closing the clamping lever, make sure that the desired clamping position number (1, 2 or 3) is visible through the borehole of the chisel, and that the clamping lever is completely closed. (Fig. 6)





Insert the chisels into the designated opening of the guide slots of the Chuck by hand until bone contact is made (Fig. 7).

Fig. 7



Remark

Make sure that the chisels are guided through both consecutively arranged guide slots of the Chuck. The chisels must not rest on the stem coating but must have direct contact with the bone before being driven in with light hammer strokes (Fig. 8).

Fig. 8



Subsequently, drive the Prestarter Chisels in along the prosthesis with light hammer strokes (Fig. 9).

Remark

Exercise caution when driving in the chisels, to minimise the risk of fissures.





Fig. 10

As soon as the Handle with the clamped chisel reaches the end face of the stem cone (or shortly before, Fig. 10), the clamping position of the chisel must be extended (clamping position 2 or 3).

Remark

Only the Starter Chisels have a mechanical stop and can be driven in down to the end position.



For re-clamping, completely open the clamping lever of the Handle, and insert the chisel in the next position. Then close the clamping lever again (Fig. 11).



The clamping lever must sit completely flat on the Handle.

Remark

To check whether the chisel runs along the stem, radiography using an image intensifier can be performed during the chiselling process.

Fig. 11



Once preparation along the lateral side of the stem has been completed with the Prestarter Chisels, continue with the Starter Chisels in the same way (Fig. 12).

Remark

Use of the medial chisels (Prestarter, Starter and Final) is only recommended in cases of good bone quality, and only when the stem cannot be loosened after preparation with the lateral, ventral and dorsal chisels.

Fig. 12



Once the Starter Chisel has reached the stop on the Chuck (Fig. 13), it is recommended to mark the bone with the axis along which this chisel was driven, for later use with the Final Chisels.

Fig. 13



Fig. 14

Remove the Starter Chisel from the bone with light backward strokes. The strikeout bolt can assist with this; screw it into the top or bottom of the Handle and strike against it with the mallet (Fig. 14).



It is recommended to start laterally with the chisels provided for this purpose and then continue along the ventral and dorsal contours in the same manner as the procedure described above (Fig 15).

Note: The ventral/dorsal chisels are used as ventral or dorsal chisels depending on whether surgery is performed on a left or right hip (see picture below and chapter 4.3 for a detailed overview).





Fig. 16

After removing the chisels, detach the Chuck from the cone with the screwdriver and remove it (Fig. 16).



Fig. 17



Insert the Final Chisel lateral into the Handle. The correct size assignment is shown in Table 1.

Make sure that the chisels are always placed in the Handle with the numbers facing upwards (the numbers must be legible when the clamping lever is closed, Fig. 17).

optimys sizes	Chisel size
1, 2, 3	Final Chisel lateral 01–03
4, 5, 6	Final Chisel lateral 04–06
7,8,9	Final Chisel lateral 07–09
10, 11, 12	Final Chisel lateral 10–12

Table 1: The size assignment is identical for all the chisels (dorsal, ventral and medial).

Remark

Each chisel can be used for three stem sizes. All the chisels have a size marking that indicates how deep the chisel may be driven into the bone for the respective stem size. The reference point is the proximal lateral corner of the stem (Fig. 18).

The Handle has three different clamping positions. The chisels are always inserted in clamping position 1 at the beginning of the preparation of each respective contour, to ensure sufficient stability of the chisels. The number <1> must be legible through the round borehole in the chisels.





When closing the clamping lever, make sure that the required clamping position (1, 2 or 3) is visible through the borehole in the chisel (Fig. 19).

Fig. 19



Fig. 20

Place the Final Chisel on the implant-bone interface, then check that its axial alignment and position are the same as the Prestarter and Starter Chisels (Fig. 20).

Remark

This is necessary to optimise use of the guiding path created with the previous chisels and successfully loosen the prosthesis-bone-interface distally.

Subsequently, drive the Final Chisel in and along the stem in axial alignment using light hammer strokes.

Remark

Exercise caution when driving in the chisels, to minimise the risk of fissures.



As soon as the Handle with the clamped chisel reaches the end face of the stem cone (or shortly before, Fig. 19), the clamping position of the chisel must be extended (clamping position 2 or 3), to prevent damage to the Handle.

For re-clamping, completely open the clamping lever of the Handle, insert the chisel in the next position, then close the clamping lever again.



The clamping lever must sit completely flat on the Handle.

Remark

To check whether the chisel runs along the stem, radiography using an image intensifier can be performed during the chiselling process.

Fig. 21



Once the chisel has been driven in down to the marked stopping point (Fig 22), remove the chisel from the bone with light backward hammer strokes, and detach the Handle from the chisel.





Other subsequent Final Chisels are driven in the same way.

It is recommended to start laterally with the chisels provided for this purpose and then work on the ventral and dorsal contours. (Fig 23)

Remark

The ventral/dorsal chisels are used as ventral or dorsal chisels depending on whether surgery is performed on a left or right hip (see chapter 4.3 for a detailed overview).

Fig. 23

After each of the required Final Chisels has been used and the stem has been sufficiently loosened, a universal extraction instrument can be attached to the prosthesis neck, and the stem may be removed from the femur with strong hammer strokes.

After removal of the prosthesis, check the femur for possible fissures, fix it if necessary, and remove foreign debris, before implanting a new prosthetic system using the implant-specific surgical technique.

4. Instruments



4.1 Basic instr. set for optimys stem removal 51.34.1173A – Configuration

ltem no.	Description
400.001	Handle
ltem no.	Description
411.921	Chuck standard 01
411.922	Chuck standard 02
411.923	Chuck standard 03
411.924	Chuck standard 04
411.925	Chuck standard 05
411.926	Chuck standard 06
411.927	Chuck standard 07
411.928	Chuck standard 08
411.929	Chuck standard 09
411.930	Chuck standard 10
411.931	Chuck standard 11
411.932	Chuck standard 12
411.901	Chuck lateral 01
411.902	Chuck lateral 02
411.903	Chuck lateral 03
411.904	Chuck lateral 04
411.905	Chuck lateral 05
411.906	Chuck lateral 06
411.907	Chuck lateral 07
411.908	Chuck lateral 08
411.909	Chuck lateral 09
411.910	Chuck lateral 10
411.911	Chuck lateral 11
411.912	Chuck lateral 12
Item no.	Description

item no.	Description
400.101	Fixation screw M4x11.5

4.2 Detail Instrument List

(re-usable instruments)







ltem no.	Description
400.001	Handle

Itom no	Description
1141 024	Description
411.921	
411.922	Chuck standard 02
411.923	Chuck standard 03
411.924	Chuck standard 04
411.925	Chuck standard 05
411.926	Chuck standard 06
411.927	Chuck standard 07
411.928	Chuck standard 08
411.929	Chuck standard 09
411.930	Chuck standard 10
411.931	Chuck standard 11
411.932	Chuck standard 12
411.901	Chuck lateral 01
411.902	Chuck lateral 02
411.903	Chuck lateral 03
411.904	Chuck lateral 04
411.905	Chuck lateral 05
411.906	Chuck lateral 06
411.907	Chuck lateral 07
411.908	Chuck lateral 08
411.909	Chuck lateral 09
411.910	Chuck lateral 10
411.911	Chuck lateral 11
411.912	Chuck lateral 12
Iteration and a	Description (

Item no.	Description
400.101	Fixation screw M4x11.5

4.3 Detail Instrument List

(single use instruments)



Item no.	Description
411.2015	Prestarter Chisel lateral
411.202S*	Prestarter Chisel medial
المحمد والمكرم والع	I Development and the later of the second

*Use of the medial Prestarter chisel is only recommended in cases of good bone quality.



ltem no.	Description
411.210S	Starter Chisel lateral 01–03
411.2115	Starter Chisel lateral 04–06
411.2125	Starter Chisel lateral 07-09
411.2135	Starter Chisel lateral 10–12
411.214S*	Starter Chisel medial 01–03
411.215S*	Starter Chisel medial 04–06
411.216S*	Starter Chisel medial 07-09
411.217S*	Starter Chisel medial 10-12
411.2145* 411.2155* 411.2165* 411.2175*	Starter Chisel medial 01–03 Starter Chisel medial 04–06 Starter Chisel medial 07–09 Starter Chisel medial 10–12

*Use of the medial Starter Chisel is only recommended in cases of good bone quality.

ltem no.	Description	right hip		left hip	
		ventral	dorsal	ventral	dorsal
411.4015	Starter Chisel ventral/dorsal 01-03	Х	_	-	Х
411.4025	Starter Chisel ventral/dorsal 01-03	-	Х	Х	-
411.403S	Starter Chisel ventral/dorsal 04-06	Х	_	-	Х
411.404S	Starter Chisel ventral/dorsal 04–06	-	Х	Х	-
411.405S	Starter Chisel ventral/dorsal 07-09	Х	_	-	Х
411.406S	Starter Chisel ventral/dorsal 07-09	-	Х	Х	-
411.407S	Starter Chisel ventral/dorsal 10-12	Х	_	-	Х
411.408S	Starter Chisel ventral/dorsal 10-12	-	Х	Х	-

ltem no.	Description
411.2185	Final Chisel lateral 01–03
411.2195	Final Chisel lateral 04–06
411.220S	Final Chisel lateral 07–09
411.2215	Final Chisel lateral 10–12
411.222S*	Final Chisel medial 01–03
411.223S*	Final Chisel medial 04–06
411.224S*	Final Chisel medial 07–09
411.225S*	Final Chisel medial 10–12

*Use of the medial final Chisels is only recommended in cases of good bone quality.

ltem no.	Description	right hip		left hip	
		ventral	dorsal	ventral	dorsal
411.4095	Final Chisel ventral/dorsal 01–03	Х	-	-	Х
411.4105	Final Chisel ventral/dorsal 01–03	-	Х	Х	-
411.4115	Final Chisel ventral/dorsal 04–06	Х	-	-	Х
411.4125	Final Chisel ventral/dorsal 04–06	-	Х	Х	-
411.4135	Final Chisel ventral/dorsal 07–09	Х	-	-	Х
411.414S	Final Chisel ventral/dorsal 07–09	-	Х	Х	-
411.415S	Final Chisel ventral/dorsal 10-12	Х	-	_	Х
411.4165	Final Chisel ventral/dorsal 10–12	-	Х	Х	-







5. Symbols



Notes

Notes

GOMINA AG

Distributed by



Australia	Mathys Orthopaedics Pty Ltd Artarmon, NSW 2064 Tel: +61 2 9417 9200 info.au@mathysmedical.com	Italy	Mathys Ortopedia S.r.l. 20141 Milan Tel: +39 02 4959 8085 info.it@mathysmedical.com
Austria	Mathys Orthopädie GmbH 2351 Wiener Neudorf Tel: +43 2236 860 999 info.at@mathysmedical.com	Japan	Mathys KK Tokyo 108-0075 Tel: +81 3 3474 6900 info.jp@mathysmedical.com
Belgium	Mathys Orthopaedics Belux N.VS.A. 3001 Leuven Tel: +32 16 38 81 20 info.be@mathysmedical.com	New Zealand	Mathys Ltd. Auckland Tel: +64 9 478 39 00 info.nz@mathysmedical.com
France	Mathys Orthopédie S.A.S 63360 Gerzat Tel: +33 4 73 23 95 95 info.fr@mathysmedical.com	Netherlands	Mathys Orthopaedics B.V. 3001 Leuven Tel: +31 88 1300 500 info.nl@mathysmedical.com
Germany	Mathys Orthopädie GmbH «Centre of Excellence Sales» Bochum 44809 Bochum Tel: +49 234 588 59 0 sales.de@mathysmedical.com	P. R. China	Mathys (Shanghai) Medical Device Trading Co., Ltc Shanghai, 200041 Tel: +86 21 6170 2655 info.cn@mathysmedical.com
	«Centre of Excellence Ceramics» Mörsdorf 07646 Mörsdorf/Thür. Tel: +49 364 284 94 0 info.de@mathysmedical.com	Switzerland	Mathys (Schweiz) GmbH 2544 Bettlach Tel: +41 32 644 1 458 info@mathysmedical.com
	«Centre of Excellence Production» Hermsdorf 07629 Hermsdorf Tel: +49 364 284 94 110 info.de@mathysmedical.com	United Kingdom	Mathys Orthopaedics Ltd Alton, Hampshire GU34 2QL Tel: +44 8450 580 938 info.uk@mathysmedical.com

Local Marketing Partners in over 30 countries worldwide...