


vitamys[®]

The E-Factor Makes the Difference



vitamys[®]

Introduction

vitamys is the highly cross-linked polyethylene (HXLPE) from Mathys. Designed for durability and longevity, vitamys meets the highest expectations of a modern orthopaedic implant material. Enriched with vitamin E, this HXLPE heralds the next generation of HXLPE.

Material benefits

- **High oxidation resistance**
- **High ageing resistance**
- **High wear resistance**
- **Excellent mechanical properties**



Radiation cross linking and thermal treatment of ultra-high molecular weight polyethylene (UHMWPE) have aroused intense scientific and commercial interest within the orthopaedic community since the late 1990s. Since 1998, highly cross-linked polyethylenes (HXLPE) are in clinical use. A reduction in wear of 25 % to 83 % after 5 years of implantation has been reported in literature¹. Despite the improvement over conventional polyethylene, previous generations of HXLPE are characterized by some major disadvantages:

- High-energy radiation destroys molecular chains, giving rise to unstable, chemically active areas of damage in the molecules, known as free radicals.
- Subsequent heat treatment causes these radicals to «pair up» more easily with other molecule chains and react with them.
- Heat treatment leads to optimized oxidation resistance but can weaken the mechanical properties of the polyethylene.

The trade-off between high oxidation resistance and high mechanical properties has now been solved with vitamys®.

vitamys®

A European development

Patented by Mathys – vitamin E additive replaces heat treatment

The patent-protected procedure by Mathys eliminates the need for heat treatment to quench the free radicals.

In this procedure, the oxidation resistance of HXLPE is assured by the inclusion of the well-known antioxidant vitamin E.

Vitamin E

Vitamin E is a collective term for the group of tocopherols, of which alpha-tocopherol (Fig. 1) shows the best properties as an antioxidant². It exists as a natural substance in the human body and also in dairy products; for example in nuts, oils and some fruits such as avocados. The vitamin E content of vitamys is below the recommended daily amount of 10–15 mg³, max. 400 mg. For example, 200g almonds contain more vitamin E than a vitamys implant. Therefore no systemic reactions are expected.

The highly cross-linked polyethylene of Mathys – vitamys® – is manufactured from UHMWPE Chirulen® 1020 and endowed with 0.1 % of synthetic vitamin E. It is homogeneously mixed with vitamin E, sintered under pressure, and highly cross-linked.




Fig. 1 Molecular structure of α -tocopherol

Homogeneous penetration of vitamin E in the polyethylene matrix

Vitamin E is fully integrated in the polyethylene, thanks to homogeneous mixing with UHMWPE Chirulen 1020 at a very early stage in the production process.⁴

The advantages

- Excellent protection against oxidation already established in the raw material
- No separate heat treatment required to diffuse vitamin E into the polyethylene
- No risk of uneven concentration profiles of vitamin E in the polyethylene⁵

Vitamin E is a stabilizer to prevent the polyethylene from oxidation – the smart solution!

A European development

Function of vitamin E during the cross linking process

Optimal thermal treatment and the addition of vitamin E make vitamys one of the most durable HXLPE available today. The addition of vitamin E as anti-oxidant is a method to preserve mechanical properties and ageing resistance. During and after the cross linking process vitamin E binds oxygen and allows the free radicals to reconnect the carbon-hydrogen (CH_2) chains (Fig. 2).

Unlike most HXLPE on the market, Mathys uses only a stress relieving thermal treatment far below melting temperature to ensure the dimensional stability of the material⁶.





Fig. 2

With vitamys[®], no extensive heat treatment is required to preserve high mechanical properties.



High oxidation, ageing and wear resistance

High oxidation and ageing resistance

Accelerated ageing through increased temperatures and oxygen content is used to simulate oxidation and ageing over many years of shelf storage. For example, 15 days accelerated ageing under O₂ (70°C, 5 bar) correlates to 10 years in vivo.⁷

In the in-vitro testing as shown in Fig. 3 below, vitamys proved to be long-term oxidation and ageing resistant.⁸ With the addition of vitamin E, vitamys demonstrated oxidation resistance for a period of 40 years or more under physiological conditions (extrapolation).




Fig. 3 Accelerated ageing test of vitamys under O₂ (5 bar at 70°C tested at 0, 14, 24, 42 days)

The comparison of conventional UHMWPE made from GUR 1020 and vitamys demonstrated the effect of vitamin E:⁸

- vitamys shows no detectable oxidation tested under O₂ at any time
- vitamys shows 10 times lower oxidation index compared to GUR 1020 tested under H₂O₂

High wear resistance

In hip simulator tests (according to ISO 14242-1:2012) vitamys proves to significantly reduce wear compared to UHMWPE^{9,10}. Wear rate of vitamys remained at constant low level even using different head materials and diameters⁹.




Fig. 4

Wear reduction with vitamys *in vivo*

After 1–2 years the RM Pressfit vitamys Cup clinically showed a significantly lower wear rate than standard UHMWPE¹¹ and confirms the positive results seen in the simulator studies.




Fig. 5

Excellent mechanical properties

Excellent mechanical properties

vitamys is a highly cross-linked polyethylene used in total hip replacement that is equivalent to the best PE type 1 classification (according to ISO), meeting all requirements for yield strength, ultimate tensile strength and elongation at break. These mechanical properties are most important for long-term performance of the material.

The figure below shows vitamys in comparison to the different PE types according to ISO. vitamys shows excellent values. Other HXLPE may only meet PE Type 2 or 3 categories due to decreased tensile strength or elongation.

Although second generation HXLPE materials are convincing in terms of wear resistance, the mechanical properties of other HXLPE have not yet met the most stringent requirements. Extensive heat treatment which is necessary to get rid of free radicals weakens the mechanical properties of polyethylene¹².




Fig. 6 vitamys in comparison to ISO standards¹³

Solutions with vitamys®

vitamys is a high wear resistant alternative to hard/hard bearings for young and active patients. Used in combination with our ceramic heads there is no risk of increased ion levels or risk of chipping/fracture of the inlay – a truly «forgiving» combination.

RM Pressfit vitamys® Cup

The RM Pressfit vitamys implant is a cementless, 100 % PE acetabular cup and is the first cup of this type to combine the advantages of low rigidity (similar to cancellous bone, thus preventing stress shielding) with resistance to wear and ageing, and clinically proven biologic anchorage.


The benefits of vitamys allow a reduction in the polyethylene wall thickness and larger articulation diameters which can be used with smaller cups: a 32 mm articulation with a cup diameter of 48 mm and a 36 mm articulation with 52 mm cups.

The RM Pressfit vitamys Cup is characterised by

- Extraordinary combination of implant design and material properties
- High wear and ageing resistance
- Proven concept of elasticity
- Easy surgical technique, MIS adapted instruments

seleXys® vitamys® Inlays

The successful modular cup system seleXys is offered with vitamys inlays in addition to conventional UHMWPE inlays. Articulations of 28, 32 and 36 mm are available.



RM Pressfit vitamys®



seleXys® PC vitamys®

vitamys® – The E-Factor Makes the Difference

Literature

- ¹ Dorr et al, JBJS 87A, p. 1816-21, 2005; D'Antonio et al, CORR 441, p. 143-50, 2005; Triclot et al, JBJS 89B, p. 1439-45, 2007; Bragdon et al, CORR 465, p. 122-7, 2007; Röhl et al, Acta Orthop. 78, p. 739-45, 2007; Digas et al, Acta Orthop. 78, p. 746-54, 2007; Olyslaegers et al, J. Arthroplasty 23, p. 489-94, 2008; Garcia-Rey et al, JBJS 90B, p. 149-53, 2008
- ² Oral E et al, The effect of alpha-tocopherol on the oxidation and free radical decay in irradiated UHMWPE. Biomaterials 2006; 27: 5580–87.
- ³ www.diabetes-news.de/info/ernaehrung/vitamine/vitamin-e.htm
- ⁴ Patent No. WO0049079; Data on file
- ⁵ Oral E et al, Wear Resistance and Mechanical Properties of Highly Cross-linked, Ultrahigh-Molecular Weight Polyethylene Doped With Vitamin E, Journal of Arthroplasty Vol.21 No.4 2006: 580-91
- ⁶ Annealing: Annealing is a heat treatment process below melting temperature wherein a material is altered, causing changes in its properties such as strength and hardness. It is used to induce ductility, relieve internal stresses, refine the structure by making it homogeneous
- ⁷ Zurbrügg D et al, 18th European Conference on Biomaterials, 2003, Stuttgart, Germany
- ⁸ Lerf R et al, Use of vitamin E to protect cross-linked UHMWPE from oxidation, Biomaterials 2010;31: 3643–48 or doi: 10.1016/j.biomaterials.2010.01.076
- ⁹ Beck M, Lerf R, Becker R et al (2012): Oxidation prevention with vitamin E in a HXLPE isoelastic monoblock pressfit cup: Preliminary results in Knahr K (Ed.), Total Hip Arthroplasty, Springer Press, 2012)
- ¹⁰ Delfosse D et al, What Happens to the Vitamin E in a Vitamin-Stabilised HXLPE, Tribology in Total Hip and Knee Arthroplasty, Springer 2014, 197-205
- ¹¹ Highly Cross-linked Polyethylene, 14th EFORT Congress 5–8 June 2013, Istanbul
- ¹² Oral E et al, Mechanisms of decrease in fatigue crack propagation resistance in irradiated and melted UHMWPE, Biomaterials 2006;27: 917–25.
- ¹³ ISO Standard 3834-2

Australia	Mathys Orthopaedics Pty Ltd Lane Cove West, NSW 2066 Tel: +61 2 9417 9200 info.au@mathysmedical.com	Japan	Mathys KK Tokyo 108-0075 Tel: +81 3 3474 6900 info.jp@mathysmedical.com
Austria	Mathys Orthopädie GmbH 2351 Wiener Neudorf Tel: +43 2236 860 999 info.at@mathysmedical.com	New Zealand	Mathys Ltd. Auckland Tel: +64 9 478 39 00 info.nz@mathysmedical.com
Belgium	Mathys Orthopaedics Belux N.V.-S.A. 3001 Leuven Tel: +32 16 38 81 20 info.be@mathysmedical.com	Netherlands	Mathys Orthopaedics B.V. 3905 PH Veenendaal Tel: +31 318 531 950 info.nl@mathysmedical.com
France	Mathys Orthopédie S.A.S 63360 Gerzat Tel: +33 4 73 23 95 95 info.fr@mathysmedical.com	P. R. China	Mathys (Shanghai) Medical Device Trading Co., Ltd Shanghai, 200041 Tel: +86 21 6170 2655 info.cn@mathysmedical.com
Germany	Mathys Orthopädie GmbH «Centre of Excellence Sales» Bochum 44791 Bochum Tel: +49 234 588 59 0 sales.de@mathysmedical.com Hotline: +49 1801 628497 (MATHYS) «Centre of Excellence Ceramics» Mörsdorf 07646 Mörsdorf/Thür. Tel: +49 364 284 94 0 info.de@mathysmedical.com «Centre of Excellence Production» Hermsdorf 07629 Hermsdorf Tel: +49 364 284 94 110 info.de@mathysmedical.com	Switzerland	Mathys Ltd Bettlach 2544 Bettlach Tel: +41 32 644 1 644 info@mathysmedical.com
		United Kingdom	Mathys Orthopaedics Ltd Alton, Hampshire GU34 2QL Tel: +44 8450 580 938 info.uk@mathysmedical.com

Local Marketing Partners in over 30 countries worldwide ...