




For healthcare professional use only. The illustrated image does not represent a connection between the use of the medical device described, nor its performance.

Preservation in motion



Ceros TCP Putty
Kneadable bone substitute

Product information


Our Ceros TCP Putty – Features & benefits

Chemical composition

Ceros TCP Putty is a sterile, osteoconductive, resorbable bone substitute offering kneadable properties. ^{1, 2, 3}

Ceros TCP Putty is made of synthetic, monophasic beta-tricalcium phosphate (β -TCP, beta-TCP, $[\text{Ca}_3(\text{PO}_4)_2]$) granules in a fermented, non-animal-derived sodium hyaluronate carrier in granular form. Ceros TCP Putty is delivered as dry powder mixture in a cup. The Ceros TCP Putty dry powder mixture consists of 92.5 % beta-TCP granules and 7.5 % sodium hyaluronate (dry w/w-%).

A kneadable Ceros TCP Putty paste is obtained by adding and mixing with a liquid component intra-operatively. The liquid component, which can consist of a sterile physiological salt solution, autologous blood or bone marrow or blood- or bone-marrow-derived cellular material, is not provided, but must be procured intra-operatively. The amount of liquid to be added to the respective powder mixture is listed in the instructions for use (Table 1, chapter 9.2) and on the product label.



Sterile Ceros TCP Putty in dry state: Ceros granules (β -TCP) and hyaluronic acid (HyA)

Synthetic β -tricalcium phosphate (ASTM F 1088)

- Mathys Ltd Bettlach has more than 35 years of clinical experience with synthetic β -TCP
- Synthetic, biocompatible, osteoconductive and resorbable bone graft substitute^{1, 2, 3}
- No risk of disease transmission

Fermented hyaluronic acid (EN ISO 22442 1-3)

- Non-animal derived hyaluronic acid in granular form, obtained by fermentation
- No risk of disease transmission
- Hyaluronic acid is naturally occurring in the body and one of the main components of the extracellular matrix

Remodeling into bone

The monophasic beta-tricalcium phosphate granules have defined porosity and are remodeled by cellular activity and as a rule replaced by autologous bone in vivo within 6–18 months. The sodium hyaluronate is enzymatically metabolized in vivo. Depending on the patient's constitution (sex, metabolism, smoking status) and age, and the location and size of the bone defect or intended bony fusion, the remodeling time may vary.^{1, 3}


Adding autologous blood or bone marrow or blood- or bone-marrow-derived cellular material or morselized bone to synthetic bone substitutes is state of the art to enhance their osteointegration.^{4, 5, 6}

Performance and indication

Intended use and indications

Ceros TCP Putty is used in bone defect filling in periodontal, oral and maxillofacial surgery, including filling of extraction socket. It also is used to increase bone volume with bone augmentation on the sinus floor and/or alveolar ridges using the guided bone regeneration technique.

For detailed description of preparation and indications/contraindications please read the instructions for use.



Evidence from clinical studies

No.	Field of application	Product	Reference
1	Sinus Floor Augmentation	Ceros TCP Putty	Foitzik et al (2015) ¹ Kluk (2013) ² Knabe et al. (2017) ⁷
2	Horizontal Augmentation, dehiscence, in combination with GBR	Ceros TCP granules* with autologous bone	Merli 2015_a ⁸
	Periodontal defects in combination with GBR	Ceros TCP granules* with autogenous platelet concentrate	Moder 2012 ⁹
	Lateral mandibular defect	Ceros TCP granules* with autologous bone	Merli 2013 ¹⁰
3	Alveolar defect	Ceros TCP granules*, 0.7–1.4 mm	Markwalder 1987 ¹¹

*clinical data based on equivalence product Ceros TCP Granules

Evidence from animal studies

No.	Field of application	Product	Reference
4	Mandibular bone defect in minipigs	Ceros TCP granules*, 0.7–1.4 mm	Buser 1998 ¹² Jensen 2005 ¹³
	Mandibular bone defect in minipigs	Ceros TCP granules*, 0.7–1.4 mm with blood or platelet concentrate	Jensen 2006 ¹⁴
	Mandibular bone defect in sheep	Ceros 82 granules*, 0.7–1.4 mm	Gatti 1990 ¹⁵
	Mandibular ridge augmentation in combination with GBR in dog	Ceros TCP granules*, 0.7–1.4 mm	Von Arx 2001 ¹⁶


*clinical data based on equivalence product Ceros TCP Granules

Clinical application and order information

Clinical application

- 1) Open the cup containing the Ceros TCP Putty dry powder mixture.
- 2) Fill in complete amount of liquid into blister (**Fig. A**).
- 3) Mix liquid with powder during approximately 1 minute (**Fig. B**).
- 4) Wait after mixing for at least 1 minute. Do not wait longer than for 2 hours.
- 5) Take out Ceros TCP Putty using sterile gloves or using a sterile spatula or other sterile accessories (**Fig. C**).
- 6) Knead and shape Ceros TCP Putty between your fingers wearing sterile gloves.
- 7) Shape Ceros TCP Putty (**Fig. D**).
- 8) Apply pre-shaped Ceros TCP Putty onto the desired location.
- 9) Close bone defect using appropriate procedures.
- 10) No additional waiting time is required, as Ceros TCP Putty does not set in situ after implantation (no cement reaction).

For detailed information on the clinical application including the amount of liquid to be added to the respective powder mixture please read Instruction for Use.



Order information

Ceros TCP Putty

Item no.	Description
42.34.0101	Ceros TCP Putty, 0.5 cc
42.34.0102	Ceros TCP Putty, 1.0 cc
42.34.0103	Ceros TCP Putty, 2.5 cc
42.34.0104	Ceros TCP Putty, 5.0 cc

Material: β -tricalcium phosphate, non-animal derived sodium hyaluronate

To learn more about Ceros TCP Putty or the Ceros portfolio in general please contact your local Mathys representative directly.

References

- ¹ Foitzik J., Mensing M., Wimmers S., Lates M., et al. Sinusbodenelevation mit einer optimierten Zubereitung von phasenreinem Beta-Tricalciumphosphat – eine Multicenterstudie. *Dent Implantol*, 2015. 19(3): p. 168-171.
- ² Kluk E. Einfluss eines Trikalziumphosphatgranulats mit einheitlicher Korngröße sowie eines Trikalziumphosphatputtys mit einem Gemisch von Granula unterschiedlicher Korngröße auf die Osteogenese und Volumenstabilität nach beidseitiger Sinusbodenaugmentation – ein Vergleich anhand einer randomisierten klinischen Studie mit Split-Mouth-Design. Dissertation, 2013.
- ³ Stiller M., Kluk E., Bohner M., Lopez-Heredia M.A., et al. Performance of β -tricalcium phosphate granules and putty, bone grafting materials after bilateral sinus floor augmentation in humans. *Biomaterials*, 2014. 35(10): p. 3154-3163.
- ⁴ Becker S.M., O.;Ponomarev, I.;Stoll, T.;Rahn, B.;Wilke, I. Osteopromotion by a beta-tricalcium phosphate/ bone marrow hybrid implant for use in spine surgery. *Spine (Phila Pa 1976)*, 2006. 31(1): p. 11-7.
- ⁵ Stoll T., Maissen O., Meury T., Becker S. New Aspects in Osteoinduction. *Materialwissenschaft und Werkstofftechnik*, 2004. 35(4): p. 198-202.
- ⁶ Muschik M.L., R.;Halbhubner, S.;Bursche, K.;Stoll, T. Beta-tricalcium phosphate as a bone substitute for dorsal spinal fusion in adolescent idiopathic scoliosis: preliminary results of a prospective clinical study. *Eur Spine J*, 2001. 10 Suppl 2: p. S178-84.
- ⁷ Knabe C., Adel Khattab D., Kluk E., Struck R., et al. Effect of a Particulate and a Putty-Like Tricalcium Phosphate-Based Bone-grafting Material on Bone Formation, Volume Stability and Osteogenic Marker Expression after Bilateral Sinus Floor Augmentation in Humans. *Journal of Functional Biomaterials*, 2017. 8(3): p. 31.
- ⁸ Merli M., Moscatelli M., Mariotti G., Pagliaro U., et al. Membranes and Bone Substitutes in a One-Stage Procedure for Horizontal Bone Augmentation: A Histologic Double-Blind Parallel Randomized Controlled Trial. *Int J Periodontics Restorative Dent*, 2015. 35(4): p. 271-81 / 463-71.
- ⁹ Moder D., Taubenhansl F., Hiller K.A., Schmalz G., et al. Influence of autogenous platelet concentrate on combined GTR/graft therapy in intrabony defects: a 7-year follow-up of a randomized prospective clinical split-mouth study. *J Clin Periodontol*, 2012. 39(5): p. 457-65.
- ¹⁰ Merli M., Moscatelli M., Mazzoni A., Merli M., et al. Lateral bone augmentation applying different biomaterials: A clinical and histological evaluation of a case report. *Zeitschrift für Zahnärztliche Implantologie*, 2013(1): p. 70-79.
- ¹¹ Markwalder TM. Füllung von Alveolen nach Zahnextraktionen mit Ceros TCP-Granulat. *Case Report* 1987
- ¹² Buser D. et al., Evaluation of filling materials in membrane-protected bone defects. A comparative histomorphometric study in the mandible of miniature pigs. *Clin. Oral Implants Res*. 1998 Jun; 9(3): 137-40
- ¹³ Jensen SS, Broggini N, Weibrich G, Hjorting-Hansen E, Schenk R, Buser D. Bone regeneration in standardized bone defects with autografts or bone substitutes in combination with platelet concentrate: a histologic and histomorphometric study in the mandibles of minipigs. *Int J Oral Maxillofac Implants* 2005;20:703-12.
- ¹⁴ Jensen SS, Broggini N, Hjorting-Hansen E, Schenk R, Buser D. Bone healing and graft resorption of autograft, anorganic bovine bone and beta-tricalcium phosphate. A histologic and histomorphometric study in the mandibles of minipigs. *Clin Oral Implants Res* 2006;17:237-43.
- ¹⁵ Gatti AM, Zaffe D, Poli GP. Behaviour of tricalcium phosphate and hydroxyapatite granules in sheep bone defects. *Biomaterials* 1990;11:513-7.
- ¹⁶ von Arx T, Cochran DL, Hermann JS, Schenk RK, Higginbottom FL, Buser D. Lateral ridge augmentation and implant placement: an experimental study evaluating implant osseointegration in different augmentation materials in the canine mandible. *Int J Oral Maxillofac Implants* 2001;16: 343-54.

Australia	Mathys Orthopaedics Pty Ltd Lane Cove West, NSW 2066 Tel: +61 2 9417 9200 info.au@mathysmedical.com	Italy	Mathys Ortopedia S.r.l. 20141 Milan Tel: +39 02 5354 2305 info.it@mathysmedical.com
Austria	Mathys Orthopädie GmbH 2351 Wiener Neudorf Tel: +43 2236 860 999 info.at@mathysmedical.com	Japan	Mathys KK Tokyo 108-0075 Tel: +81 3 3474 6900 info.jp@mathysmedical.com
Belgium	Mathys Orthopaedics Belux N.V.-S.A. 3001 Leuven Tel: +32 16 38 81 20 info.be@mathysmedical.com	New Zealand	Mathys Ltd. Auckland Tel: +64 9 478 39 00 info.nz@mathysmedical.com
France	Mathys Orthopédie S.A.S 63360 Gerzat Tel: +33 4 73 23 95 95 info.fr@mathysmedical.com	Netherlands	Mathys Orthopaedics B.V. 3001 Leuven Tel: +31 88 1300 500 info.nl@mathysmedical.com
Germany	Mathys Orthopädie GmbH «Centre of Excellence Sales» Bochum 44809 Bochum Tel: +49 234 588 59 0 sales.de@mathysmedical.com «Centre of Excellence Ceramics» Mörsdorf 07646 Mörsdorf/Thür. Tel: +49 364 284 94 0 info.de@mathysmedical.com «Centre of Excellence Production» Hermsdorf 07629 Hermsdorf Tel: +49 364 284 94 110 info.de@mathysmedical.com	P. R. China	Mathys (Shanghai) Medical Device Trading Co., Ltd Shanghai, 200041 Tel: +86 21 6170 2655 info.cn@mathysmedical.com
		Switzerland	Mathys (Schweiz) GmbH 2544 Bettlach Tel: +41 32 644 1 458 info@mathysmedical.com
		United Kingdom	Mathys Orthopaedics Ltd Alton, Hampshire GU34 2QL Tel: +44 8450 580 938 info.uk@mathysmedical.com

Local Marketing Partners in over 30 countries worldwide ...

Distributor

