**Processing instructions
for instruments
with flexible shafts**

Preservation in motion

For healthcare professional use only. The illustrated image does not represent a connection between the use of the medical device described, nor its performance.*

*Building on our heritage
Moving technology forward
Step by step with our clinical partners
Towards a goal of preserving mobility*

Preservation in motion

As a Swiss company, Mathys is committed to this guiding principle and pursues a product portfolio with the goal of further developing traditional philosophies with respect to materials or design in order to address existing clinical challenges. This is reflected in our imagery: traditional Swiss activities in conjunction with continuously evolving sporting equipment.

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Table of contents

1.	Scope	19
2.	Overview of the instruments with flexible shafts	19
2.1	Instruments with flexible shafts	19
3.	Terminology (EN ISO 17664)	20
3.1	Treatment chemical	20
3.2	Cleaning	20
3.3	Disinfection	20
3.4	Manual cleaning	20
3.5	Manufacturer	20
3.6	Reprocessing	20
3.7	Processor	20
3.8	Sterilisation medium	20
3.9	Sterile	20
3.10	Sterilisation	20
3.11	Validation	20
3.12	Certification	20
3.13	Washer/Disinfector	20
4.	Information to be provided by the manufacturer	20
4.1	Instructions for reprocessing	20
4.2	Limitations and restrictions of reprocessing	20
4.2.1	Warnings and precautions	20
4.2.2	Limitations	21
4.3	Preparation at the place of use – processing during and immediately after use	21
4.4	Preparation before cleaning	22
4.5	Cleaning	24
4.5.1	Instructions for manual pre-cleaning	24
4.5.2	Instructions for mechanical cleaning and disinfection	24
4.6	Maintenance/care	25
4.7	Function testing	25
4.8	Packaging	28
4.9	Sterilisation	28
4.10	Storage	29
4.11	Efficacy of the treatment process	29
4.12	Service life/Number of processing cycles	30
5.	Symbol	30
5.1	Symbol according to Mathys Ltd Bettlach	30
6.	Customer information	30
7.	Annex – quick overview	31
7.1	Manual pre-cleaning	31
7.2	Mechanical cleaning (WD)	31
7.3	Sterilisation with saturated steam	31

Processing instructions for instruments with flexible shafts

1. Scope

These Processing Instructions in accordance with the requirements of EN ISO 17664 apply to instruments that are repeatedly used and therefore reprocessed, as well as to medical devices that are sold non-sterile, but to be used sterile. Flexible shaft for drill bits (resterilisable surgical instruments) by Mathys Ltd Bettlach belong to this group.

2. Overview of the instruments with flexible shafts

2.1 Instruments with flexible shafts

Item no.

3.14.251 RM Classic flex. countersink, 1st gen.

Item no.

3.14.256 RM Classic flex. countersink 2nd gen.

Item no.

55.02.1901 RM Classic counters. reamer flex. 3rd gen.

Item no.

3.14.252 RM Classic flex. shaft, 1st gen.

Item no.

3.14.257 RM Classic flex. ream. shaft 46–52 2 gen.

Item no.

3.14.258 RM Classic flex.ream.shaft 54–64 2 gen.

Item no.

3.14.259 RM Classic flex. ream. shaft 60–68 2 gen.

Item no.

55.02.1903 RM Classic reaming shaft, flex. 3rd gen.

Item no.

3.40.501 RM Classic flexible shaft 1st gen.

Item no.

3.14.545 Shaft, flexible

Item no.

56.02.6184– Medullary reamer flex. 8–18

56.02.6194

Item no.

5502.00.2 seleXys flexible shaft

Item no.

502.08.02.00.0 Affinis glenoid drill bit

Item no.

51.34.0940 Flexible shaft for drill bits

A assembled; **B** disassembled

3. Terminology (EN ISO 17664)

3.1 Treatment chemical

Formulation of compounds intended for use in reprocessing.

3.2 Cleaning

Removal of contamination from an object to an extent that is necessary for further processing or for the intended use.

3.3 Disinfection

Method used to reduce the number of viable microorganisms on a product to a predetermined level that is suitable for the further handling or use of the product.

3.4 Manual cleaning

Cleaning without use of a washer/disinfector.

3.5 Manufacturer

Organisation that is responsible for the design, production, packaging and labelling of a medical device before it is marketed under its own name, regardless of whether these operations are carried out by this organisation itself or by a third party commissioned accordingly.

3.6 Reprocessing

Any activity, including cleaning, disinfection and sterilisation, that is necessary to prepare a new or used medical device for its intended use.

3.7 Processor

Organisation responsible for carrying out the actions necessary to prepare a new or used medical device for its intended use.

3.8 Sterilisation medium

Treatment chemical capable of destroying microorganisms, including viruses, when used at the correct dilution/dosage and applied for the recommended contact time.

3.9 Sterile

Free of any viable microorganisms.

3.10 Sterilisation

Validated method to kill all types of viable microorganisms on the product.

3.11 Validation

Documented procedure for collecting, recording and interpreting the results needed to furnish evidence that a

process consistently yields products that meet the specifications.

3.12 Certification

Confirmation by testing and furnishing of objective evidence that the specified requirements are met.

3.13 Washer/Disinfector

Device for cleaning and disinfection of medical devices and other items used in the medical, dental, pharmaceutical and veterinary field.

4. Information to be provided by the manufacturer

4.1 Instructions for reprocessing

These reprocessing instructions for instruments with flexible shafts are based on a validated manual/machine-based reprocessing method. A purely manual or purely machine-based treatment method has not been validated by Mathys Ltd Bettlach and does not lead to sufficient cleaning success.

4.2 Limitations and restrictions of reprocessing

4.2.1 Warnings and precautions

Staff who come into contact with actually or potentially contaminated surgical instruments must take the generally accepted precautionary measures regarding personal protective equipment (RKI recommendation: gown, mask, goggles or face shields, puncture-resistant gloves, shoes, shoe covers, etc.) to avoid contact with contaminated or potentially contaminated materials, instruments and products.

Particular care must be taken when handling cutting instruments (reamers, drill bits, rasp, chisels), as they pose a risk of injury to the patients on the one hand and to the staff (operating theatre and CSSD personnel) on the other.

It must be clarified in advance that the patients, as well as the staff (operating theatre and CSSD personnel), do not react with allergic reactions due to material intolerance (various steels and plastic materials) upon direct contact with instruments.

When handling instruments with sharp points or edges, special caution is advised.

For new instruments delivered to the hospital, Mathys Ltd Bettlach recommends triple cleaning before use to build up the protective oxide layer.

For the manual cleaning process, use low-foam detergents to ensure the instruments remain visible. In manual cleaning with brushes, it is recommended to keep the instruments always below the surface of the cleaning solution. This ensures that no aerosols form and splashes are avoided that might spread contaminants. To prevent accumulation of detergent residues, the detergents must be completely removed by thorough rinsing of the product surfaces.

Heavy instruments must not be placed on top of sensitive instruments.

Do not allow contaminated instruments to dry before reprocessing them. This is important because all the steps described below for cleaning and sterilisation are facilitated by preventing blood, body fluids, bones and tissue residues as well as saline or disinfectant from drying on used instruments.

The chloride and iodide ions contained in cleaning agents and disinfectants cause pitting, so the contact of the instruments with such ions must be kept as brief as possible. Then rinse thoroughly with deionised water (DI water) to remove any residues.

The condensation moisture resulting from sterilisation can be prevented by extension of the drying phase. Excessive concentrations of detergents and disinfectants, as well as strongly acidic and alkaline detergents, can damage the protective oxide layer and cause pitting.

When using such means, the concentration and exposure time recommended by the manufacturers should be strictly adhered to.

For cleaning reusable instruments, Mathys Ltd Bettlach recommends a combined manual/mechanical cleaning process with a detergent and enzymatic solution at an alkaline pH of <12.

It is extremely important to neutralise the alkaline cleaner completely and rinse it thoroughly off the instruments. In mechanical cleaning, the information provided by the manufacturers of the equipment and cleaning agents must be followed.

Only instruments by Mathys Ltd Bettlach may be used for the placement of implants by Mathys Ltd Bettlach (see the respective surgical technique); instruments by other legal manufacturers may not.

No additional lettering of any kind may be applied to the instruments.

The instruments are packaged separately and delivered in non-sterile condition. The packaging materials must be disposed of in accordance with the local and country-

specific regulations.

If these processing instructions are adhered to in terms of care and maintenance, neither damage to nor limitation of the service life of the medical devices are to be expected.

4.2.2 Limitations

Patients considered as risk patients with regard to prion diseases such as transmissible spongiform encephalopathy (TSE) and Creutzfeldt-Jakob disease (CJD/vCJD) and the associated infections must be operated on with single-use instruments.

Flexible shafts (as well as drill bits, reamers, rasps and other cutting instruments) that have been cleaned for reprocessing with alkaline detergents must be carefully examined afterwards. It must be ensured that the cutting edges are sharp for use. If this is no longer the case, please contact your local Mathys partner.

Some enzyme solutions are specifically designed for the decomposition of faeces or other organic contaminants and are not suitable for cleaning surgical instruments. For degradation of blood, body fluids and tissues, an appropriate enzyme-containing liquid cleaner should be selected.

No restrictions exist in the repeated treatment of the instruments of Mathys Ltd Bettlach, as this process has a low impact on the life of the instruments. The service life is generally defined by wear and damage from use.

Use of hard water (${}^{\circ}\text{dH}$ value >14) is to be avoided. For the initial rinsing, soft tap water is suitable. Rinse thoroughly with deionised water (DI water) to remove any residues. Ordinary water often contains high concentrations of minerals (e.g. lime or silica), which can be seen on the instrument surface as stains with sharply defined edges.

4.3 Preparation at the place of use – processing during and immediately after use

First, residual body fluids and tissues must be removed below the water surface using a specific plastic brush made of nylon. Here it should be made sure that the water used does not exceed a temperature of 45°C (113°F), since otherwise the contaminations, such as proteins in the blood and tissues, will be fixated on the instrument. Saline, blood, body fluids, tissue, bone residues or other organic particles on instruments must be removed before cleaning as soon as possible in order to avoid drying as well as corrosion.

Tip

Immersion of the used instruments into suitable enzyme solutions or into cold deionised water (DI water) facilitates cleaning. The enzyme solutions degrade the protein-containing substances, preventing blood and protein-containing materials from drying on the instruments.

Saline solutions as well as detergents and disinfectants that contain aldehydes, mercury, active chlorine, chloride, bromine, bromide, iodine or iodide are corrosive and must not be used.

Read the instructions, and contact the manufacturer in case of doubt.

Excessive concentrations of the cleaning agents, as well as strongly acidic and alkaline detergents, can damage the protective oxide layer and cause pitting. When using such agents, the concentration and exposure time recommended by the manufacturers must be strictly adhered to. Mathys Ltd Bettlach recommends use of alkaline cleaning agents with a pH < 12.

Optimal cleaning is ensured if the instruments are properly cleaned within one hour after use in the CSSD to minimise the risk of drying of substances and materials. If it is not possible to clean the instruments within this specified time, Mathys Ltd Bettlach recommends immersing the instruments into an enzyme solution or deionised water (DI water) of room temperature, or to wrap them for up to 6 hours into cloths wetted with enzyme solution or deionised water (DI water).

For prevention of contamination, the used instruments must be transported in closed or covered containers to the central sterilisation department.

Strongly contaminated solutions (blood and/or clouding) should be replaced with freshly prepared cleaning solutions.

If water without detergent additives is used for cleaning, Mathys Ltd Bettlach recommends a water temperature of no more than 45°C (113°F), since otherwise the blood will start to coagulate and hence proteins are going to be strongly fixed on the instrument, which can be removed only with major effort in mechanical cleaning.

4.4 Preparation before cleaning

In mechanical cleaning, the information provided by the manufacturers of the equipment and cleaning agents must be followed.

When using dry detergents in powder form, it must be made sure that these are completely dissolved before use to avoid discolouration or corrosion of the instruments.

Table 1: Overview of the reprocessing according to EN ISO 17664

Procedure			Reusable surgical instruments
Preparation at the place of use	Condition	Wet	<ul style="list-style-type: none"> Recommendation cold deionised water (DI water) or enzyme solution (fluid or wetted cloths)
		Drying	<ul style="list-style-type: none"> Max. 6 hours Max. 1 hour
Decontamination	Preparation	• Selection according to the cleaning and disinfection method	
	Cleaning	Manual	–
		Mechanical	–
		Ultrasound	+
		Combined manual/mechanical	+
		Alkaline (pH < 12)	+
		Acidic	–
		Neutral to slightly alkaline (pH 7 – 9.5)	–
	Rinsing	Final rinse with deionised water (DI water)	
	Disinfection ¹	Chemically max. 60°C (140°F)	–
		Thermally 90°C (194°F)	+
	Drying	T _{max} /Time	115°C (239°F)/15 min
Maintenance	Function check		Mandatory
	Care	Sterilit® I JG600	Mandatory
Sterilisation	Moist heat ²		+
	Ethylene oxide		–
	Formaldehyde		–
	Plasma		–

Legend: + validated method, – non-validated method

¹ Non-automatic disinfection process according to the «Guideline for Disinfection and Sterilisation in Healthcare Facilities, 2008»

² Preferred sterilisation method according to EN ISO 17664

In order to achieve optimum and thorough cleaning results, Mathys Ltd Bettlach prescribes a combined manual/mechanical cleaning procedure.

4.5 Cleaning

4.5.1 Instructions for Manual pre-cleaning

No.	Step	Medium
1	Manual removal of all visible dirt using a plastic brush* below the water surface with slight bending of the flexible shaft in all directions; for snap-fit/quick-release fasteners, also move them in all directions until no visible residues are present any longer.	<ul style="list-style-type: none">• Nylon brush• Tap water
2	Place a white lint-free cloth on the table and strike the flexible shaft in bent condition onto the cloth. Repeat this process in different positions. If residues form on the white cloth, repeat step 1.	<ul style="list-style-type: none">• White, lint-free cloth
3	Completely submerge the instruments into a combined detergent and enzyme solution in an ultrasonic bath having a water temperature of 35°C–45°C (95°F–113°F) and a frequency of 35–47 kHz, and sonicate them for 30 minutes.	<ul style="list-style-type: none">• Enzymatic cleaner 1 % deconex® TWIN PH10 and 1 % deconex® TWIN ZYME, (v/v) in deionised water (DI water)
4	Rinse the flexible shaft under running tap water with simultaneous slight bending of the flexible shaft in all directions and with concomitant brushing with a new plastic brush*.	<ul style="list-style-type: none">• Nylon brush• Tap water
5	Instruments with snap-fit/quick-release fasteners: Rinsing of the snap-fit/quick-release fasteners under running tap water with concomitant moving of the snap-fit/quick-release fastener in all directions and simultaneous brushing with a new plastic brush*.	<ul style="list-style-type: none">• Nylon brush• Tap water
6	Make sure that no visible residues are present any longer. For this purpose, bend the flexible shaft into various positions as well. If necessary, repeat procedure. If the outflowing water is still discoloured, repeat steps 1–5.	<ul style="list-style-type: none">• Visually

* Decontaminate brushes after use and sterilise or dispose of them. Do not use steel brushes.

4.5.2 Instructions for mechanical cleaning and disinfection

The process was carried out at Mathys Ltd Bettlach using a washer/disinfector (WD) by Miele KG (Miele Professional G7836CD) and a combined manual/machine-based cleaning process with detergent and enzyme solution deconex® TWIN PH10 and deconex® TWIN ZYME by Borer Chemie AG.

No.	Step		Medium
1	Pre-rinse	Duration: 2 minutes	• Cold tap water
2	Cleaning process	Dosage / Temperature: 0.5 % deconex® TWIN PH10 at 35°C (95°F) 0.2 % deconex® TWIN ZYME at 40°C (104°F) Duration / Temperature: 10 minutes at 55°C (131°F)	• Enzymatic cleaner 0.5 % deconex® TWIN PH10 and 0.2 % deconex® TWIN ZYME, (v/v) in deionised water (DI water)
3	Rinse I	Duration: 2 minutes Temperature: Max. 50°C (122°F)	• Tap water
4	Rinse II	Duration: 2 minutes Temperature: max. 40°C (104°F)	• Deionised water (DI water)
5	Thermal Disinfection ¹	Duration: 7 minutes Temperature: 90°C (194°F)	• Deionised water (DI water)
6	Drying	Duration: 15 minutes Temperature: 115°C (239°F)	–
7	Make sure that no visible residues are present any longer. To this end, all flexible shafts of both designs can be easily bent into various positions.		
<p>Tip <i>Bend the flexible portion of the shafts slightly and dry by striking onto a white lint-free cloth.</i></p>			

¹ Thermal disinfection according to DIN EN ISO 15883-1, Annex B, and DIN EN ISO 15883-2

4.6 Maintenance / care

After cleaning, the instruments must be completely dry and free of visible and noticeable residues. Critical areas such as handle structures, long and/or thin boreholes or blind holes, joints and complex structures must be treated with special care. To ensure that all contamination has been removed, it is of crucial importance to inspect each instrument carefully and check it for cleanliness as well as water spots (for example, lime or silicate). If any contamination should be discovered adhering to instruments, the complete manual as well as mechanical cleaning and disinfection process must be repeated immediately.

Once the instrument is visually clean, it must undergo maintenance. To this end, Mathys recommends use of a care product based on paraffin/white oil that is biocompatible, suitable for steam sterilisation, and vapour-permeable, e.g. Sterilit I JG600 (Aesculap AG). Alternative products must be oil-free, free of silicone oil-containing care products, suitable for steam sterilisation and biocompatible (see the «Red Brochure» by the AKI).

For maintenance, the instruments must be cooled to room temperature since otherwise there would be a risk of metal abrasion. The care product must be manually applied specifically, carefully, and drop by drop to hinged

or ball bearings of a snap-fit, rotating or joint mechanism and/or sliding surfaces and then distributed evenly by moving the hinges, joints, snap-fit mechanisms, or sliding surfaces. Excess care product must be removed with a lint-free cloth (the manufacturer's instructions must be observed). «Overspraying» the instruments or immersion baths is not recommended by Mathys. Plastic surfaces must not be treated with care products. Observe the expiry date indicated by the manufacturer of the care products.


4.7 Function testing

Check the instruments for damage and/or intact surfaces. Damage or wear that might impair the function of the instrument must be reported to your local Mathys Partner. The same will decide on repair and/or exchange of the instruments.

The functionality of moving parts must likewise be checked to ensure that the intended movements can be completely carried out.

Defects and their causes, as well as correct troubleshooting in case of damage, are shown in the table below.

Defect	Cause	Testing	Action
Bent shaft	<ul style="list-style-type: none"> • Improper handling • Excessive bending during application • Overloaded trays 	<ul style="list-style-type: none"> • Successful rebending possible 	<ul style="list-style-type: none"> • Continue use
		<ul style="list-style-type: none"> • Rebending impossible 	<ul style="list-style-type: none"> • Return to Mathys Ltd Bettlach or disposal, depending on the contractual relationship
Damage to the spring	<ul style="list-style-type: none"> • Improper handling • Excessive bending during application • Overloaded trays 	<ul style="list-style-type: none"> • Visual inspection detects damage to the spring 	<ul style="list-style-type: none"> • Return to Mathys Ltd Bettlach or disposal, depending on the contractual relationship
Damage/breakage of the laser-guided structure	<ul style="list-style-type: none"> • Improper handling • Excessive bending during application • Overloaded trays 	<ul style="list-style-type: none"> • Visual inspection detects damage to the laser-guided structure 	<ul style="list-style-type: none"> • Return to Mathys Ltd Bettlach or disposal, depending on the contractual relationship
Lettering no longer readable	<ul style="list-style-type: none"> • Wear through use • Fading due to cleaning process, detergent and sterilisation procedure 	<ul style="list-style-type: none"> • Item number, batch and EC marking can still be read 	<ul style="list-style-type: none"> • Continue use
		<ul style="list-style-type: none"> • Item number, batch and EC marking can no longer be read 	<ul style="list-style-type: none"> • Return to Mathys Ltd Bettlach or disposal, depending on the contractual relationship
Corrosion	<ul style="list-style-type: none"> • Unsuitable cleaning agents • Too long contact time with cleaning agents and disinfectants • Improper reprocessing 	<ul style="list-style-type: none"> • Visual examination detects corrosion 	<ul style="list-style-type: none"> • Return to Mathys Ltd Bettlach or disposal, depending on the contractual relationship

Defect	Cause	Testing	Action
Organic residues 	<ul style="list-style-type: none"> Unsuitable cleaning agents and disinfectants Too long drying of the residues before reprocessing Insufficient brushing or rinsing during the manual pre-cleaning 	<ul style="list-style-type: none"> Visual examination detects organic residues 	<ul style="list-style-type: none"> Repetition of the treatment process <p>If repeated processing is not successful:</p> <ul style="list-style-type: none"> Return to Mathys Ltd Bettlach or disposal, depending on the contractual relationship
Machine connection damaged, connector can no longer be plugged on 	<ul style="list-style-type: none"> Improper handling Overloaded trays Droppage 	<ul style="list-style-type: none"> Connection can no longer be plugged on 	<ul style="list-style-type: none"> Return to Mathys Ltd Bettlach or disposal, depending on the contractual relationship
Damage to the cutting edge 	<ul style="list-style-type: none"> Improper handling Overloaded trays Droppage 	<ul style="list-style-type: none"> Visual examination detects damage to the cutting edge 	<ul style="list-style-type: none"> Return to Mathys Ltd Bettlach or disposal, depending on the contractual relationship
Stains caused by silicate deposits or water spots caused by lime scale 	<ul style="list-style-type: none"> Rinsing with DI water not performed Carry-over of silicate-based detergents 	<ul style="list-style-type: none"> Visual examination detects stains caused by silicate deposits or water spots caused by lime scale 	<ul style="list-style-type: none"> Repetition of the treatment process
Loss of cutting power noticeable upon application	<ul style="list-style-type: none"> Regular wear and tear from use 	<ul style="list-style-type: none"> Feedback from the surgeon after use 	<ul style="list-style-type: none"> Return to Mathys Ltd Bettlach or disposal, depending on the contractual relationship

Instruments must be replaced if:

1. The surfaces look «chalky».
2. They show any signs of damage (e.g. (hairline) cracks, flaking, deformation, blistering).
3. They have excessive shape changes and/or are visibly warped.
4. The lettering, such as item no. or LOT no, is no longer legible. This likewise applies to surgical instruments that do not comprise any plastic materials and are made of steel only.

For replacement, contact your Mathys Partner.

If stains on the medical devices should be recognisable, their cause must first be ascertained. Thus, coloured spots indicate incompatibility with a process chemical or exceedance of an exposure time. White spots are often residues of lime, process chemicals or salts. Corrosion marks should not be underestimated, and affected instruments should be immediately separated from unaffected ones («flash rust» or «rust bloom»).

As damaged instruments can no longer function properly, all reprocessable instruments must be checked for proper functioning after maintenance/care but before sterilisation (see the «Red Brochure» by the AKI).

Markings on the instruments must be legible. This includes scales indicating angles, for determining implant size, of length and/or depth, and of directions such as «left» and «right». If any scales or other markings should not be legible any more, notify your local Mathys Partner promptly for assessment or for replacement of the instruments, respectively.

Please pay particular attention to the following:

1. The instrumentation must be checked for completeness.
2. The instruments in the tray must be checked for correct arrangement.
3. The instruments must be checked for damage (e.g. (hairline) cracks, deformations, changing gaps between metal and plastics, fractures, corrosion or signs of wear) and damaged surfaces. Damage or wear that might impair the function of the instrument must be reported to your local Mathys Partner. The same will decide on repair or exchange of the instruments or entire instrument trays.
4. The functionality of mobile components (e.g. hinge joints, sliding parts, moving parts, etc.) must be checked in order to ensure that the intended movement sequence can be performed and correctly.
5. Long and narrow instruments must be checked for

bending.

6. Instruments that consist of several individual components and have to be assembled for function must be checked for correct assembly and functionality after assembly.
7. Drill bits, reamers, rasps and other cutting instruments must be carefully examined for their cutting edges. It must be ensured that the cutting edges are sharp for use, and that no visible or palpable damage is present. This can be done easily using a 10–12× magnifier.
8. Instruments that are no longer functional must be returned to Mathys for repair or scrapping. Before, the instruments must undergo an entire processing cycle in order to eliminate the risk of infection.

4.8 Packaging

The packaging for the sterilisation must be suitable for the sterilisation procedure with moist heat, i.e., the permeability of the packaging for the steam must be ensured. Moreover, the packaging provides protection during transport and storage.

Therefore, Mathys Ltd Bettlach recommends double packaging of the instrument tray container.

When using sterilisation wrap, this must be free of detergent residues. Mathys Ltd Bettlach discourages use of recyclable wraps.

4.9 Sterilisation

As processor of the instruments, the hospital is responsible for performing a procedure validation of all the above steps in order to ensure complete penetration of the sterilising vapour to the instruments.

Furthermore, the hospital must take protective measures for sharp or potentially dangerous parts of instruments.

Instructions by the manufacturer of the sterilisation device must always be followed. If several sets of instruments are sterilised in one sterilisation cycle, the maximum loading of the device in accordance with the manufacturer's instructions must not be exceeded.

For optimum sterilisation, the instrument sets must be duly prepared and packaged into the appropriate trays and containers. Only in this way can the spread and penetration of steam reach all surfaces. In case of steam sterilisation, it must be ensured that the product is completely dry after sterilisation.

Steam or moist heat is the recommended sterilisation method for instruments by Mathys Ltd Bettlach (EN ISO 17664). Ethylene oxide, gas plasma and dry heat are not

recommended as sterilisation methods for the sterilisation of reusable instruments.

For sterilisation of the instruments, in any case the national recommendations and guidelines are to be complied with.

Below, the minimum sterilisation parameters are indicated that have been performed by Mathys Ltd Bettlach with a sterilisation device (Sterimed FAV6767100S or Euro-Selectomat) and validated by microbiological examination to reach a SAL (sterility assurance level) of 10^{-6} .

Sterilisation process with saturated steam

Type of cycle	Minimum temperature in °C ⁵ (°F)	Minimum sterilisation time in minutes	Minimum drying time in minutes	Minimum pressure in mbar ^{6, 7}
Pre-vacuum – pulsating vacuum (CH) ¹	134 (273)	18	20	≥3042
Pre-vacuum – pulsating vacuum (D) ²	134 (273)	5	20	≥3042
Pre-vacuum – pulsating vacuum (GB) ^{3, 4}	134 (273)	3	20	≥3042

¹ Recommendation of the World Health Organisation (WHO) for steam sterilisation of instruments with potential TSE/CJD contamination

² Hygiene requirements for the processing of medical devices, Federal Institute for Drugs and Medical Devices, 2012

³ Validated sterilisation process with a minimum sterilisation time of 3 minutes at 134°C (273°F) to achieve a Sterility Assurance Level (SAL) of 10^{-6} in accordance with EN ISO 17665-1

⁴ Validation in the original instrument case

⁵ Maximum temperature 137°C (279°F) according to SN EN 285

⁶ Pressure during the sterilisation phase according to ISO 17665-2

⁷ Pressure during the sterilisation phase at 137°C (279°F) must be ≥3318.5 mbar be according to ISO 17665-2

4.10 Storage

The packaged, sterile instruments must be stored protected from dust, insects, vermin and direct sunlight in a dry and cool room. This may be accessible only to relevant staff. The equipment for storage and transportation must be such that any disorder, overloading or falling is avoided. Sterile medical devices must never be stored directly on the floor. The instruments are used in the sequence of the receipt of the goods, and the sterile packaging of the instruments is inspected meticulously for package integrity before opening.

Each hospital must, adapted to its validated sterilisation process, specify how long the sterile-packaged instruments may be stored before the next use (DIN 58953-8).

If the packaging or a sterile wrap is torn, punctured, visibly damaged or has become wet, the instrument set must be repackaged

and sterilised.

In case of signs of open or damaged lid seals or filters on the sterilisation container, too, the instrument set must be re-sterilised and the sterile filter be replaced. For reusable filters, a careful visual inspection must be carried out.

4.11 Efficacy of the treatment process

The processing method recommended in these processing instructions has been validated. The results meet the requirements with regard to both the limit value and to the guidance value in terms of protein residues in accordance with the guideline by DGKH (Deutsche Gesellschaft für Krankenhaushygiene e. V.), DGSV (Deutsche Gesellschaft für Sterilgutversorgung e. V.) and AKI (Arbeitskreis Instrumentenaufbereitung – Working Group Instrument Reprocessing) for mechanical cleaning and thermal disinfection processes for medical devices (D 2596 F).

4.12 Service life / Number of processing cycles

The service life of the flexible shaft is generally defined by wear and damage due to use. The processor is responsible for checking for optimal functioning (e.g. cutting ability), cleanliness and defects (e.g. corrosion) prior to each use.

5. Symbol

5.1 Symbol according to Mathys Ltd Bettlach

Caution

6. Customer information

Mathys Ltd Bettlach
Robert Mathys Strasse 5
P.O. Box
2544 Bettlach
Switzerland

Tel. +41 32 644 1 644
Fax +41 32 644 1 161

info@mathysmedical.com
www.mathysmedical.com

7. Annex – quick overview

7.1 Manual pre-cleaning

No.	Step
1	Manual removal of all visible contaminations <ul style="list-style-type: none"> Below the water surface (tap water) With slight bending of the flexible shaft in all directions Using plastic brushes made of nylon Move snap-fit/quick-release fasteners into all directions as well
2	Strike the flexible shaft in bent condition (different positions) onto a white lint-free cloth . If residues form, the manual pre-cleaning must be started over again from the beginning.
3	Fill ultrasonic bath with combined detergent-enzyme solution. Completely immerse instruments and sonicate them at 35–47 kHz for 30 minutes.
4	Rinse the flexible shaft under running tap water with simultaneous slight bending of the flexible shaft in all directions and with concomitant brushing with a new plastic brush made of nylon .
5	Instruments with snap-fit/quick-release fasteners: Rinsing of the snap-fit/quick-release fasteners under running tap water with concomitant moving of the snap-fit/quick-release fastener in all directions and simultaneous brushing off with a new plastic brush made of nylon.
6	Make sure that no visible residues are present any longer. For this purpose, bend the flexible shaft into various positions as well. If necessary, repeat procedure. If the outflowing water is still discoloured, repeat steps 1–5.

7.2 Mechanical cleaning (WD)

WD connection: Place the instrument with flexible shaft into a WD cleaning tray without active water connection.

Pre-rinse	Duration: 2 minutes	• Cold tap water
Cleaning process	Duration / Temperature: 10 minutes at 55°C (131°F)	• Enzymatic cleaner 0.5 % deconex® TWIN PH10 and 0.2 % deconex® TWIN ZYME, (v/v) in deionised water (DI water)
Rinse I	Duration: 2 minutes Temperature: Max. 50°C (122°F)	• Tap water
Rinse II	Duration: 2 minutes Temperature: Max. 40°C (104°F)	• DI water
Thermal disinfection	Duration: 7 minutes Temperature: 90°C (194°F)	• DI water
Drying	Duration: 15 minutes Temperature: 115°C (239°F)	–

7.3 Sterilisation with saturated steam

Type of cycle	Minimum temperature in °C (°F)	Minimum sterilisation time in minutes	Minimum drying time in minutes	Minimum pressure in mbar
Pre-vacuum – pulsating vacuum (CH) ¹	134 (273)	18	20	≥3042
Pre-vacuum – pulsating vacuum (D)	134 (273)	5	20	≥3042
Pre-vacuum – pulsating vacuum (GB) ²	134 (273)	3	20	≥3042

¹ Recommended sterilisation process, ² validated sterilisation process

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