**MATHYS**   
a company of enovis™

**15 YEARS  
CLINICAL  
EXPERIENCE**

# Affinis Inverse

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# Literaturangaben

- <sup>1</sup> Irlenbusch U and Kohut G. Evaluation of a new baseplate in reverse total shoulder arthroplasty – comparison of biomechanical testing of stability with roentgenological follow up criteria. *Orthop Traumatol Surg Res.* 2015;101(2):185-90.
- <sup>2</sup> Kohut G, Dallmann F, Irlenbusch U. Wear-induced loss of mass in reversed total shoulder arthroplasty with conventional and inverted bearing materials. *J Biomech.* 2012;45(3):469-73.
- <sup>3</sup> Lerf R, Wuttke V, Reimelt I, Dallmann F, Delfosse D. Tribological Behaviour of the «Reverse» Inverse Shoulder Prosthesis. 7<sup>th</sup> International UHMWPE Meeting; Philadelphia 2015.
- <sup>4</sup> Irlenbusch U, Kaab M, Kohut G, Proust J, Reuther F, Joudet T. Reversed shoulder arthroplasty with inversed bearing materials: 2-year clinical and radiographic results in 101 patients. *Arch Orthop Trauma Surg.* 2015;135(2):161-9.
- <sup>5</sup> Kohut G, Reuther F, Joudet T, Kääh MJ, Irlenbusch U. Inverted-bearing reverse total shoulder arthroplasty: scapular notching does not affect clinical outcomes and complications at up to 7 years of follow-up. *J Shoulder Elbow Surg.* 2021;31(4):868-74.
- <sup>6</sup> Kääh M, Kohut G, Irlenbusch U, Joudet T, Reuther F. Reverse total shoulder arthroplasty in massive rotator cuff tears: does the Hamada classification predict clinical outcomes? *Arch Orthop Trauma Surg.* 2021;142(7):1405-11.
- <sup>7</sup> Huber J, Irlenbusch U, Kääh MJ, Reuther F, Kohut G, Judge A. Treatment effects of reverse total shoulder arthroplasty – a simple method to measure outcomes at 6, 12, 24 and 60 months for each patient. *BMC Musculoskelet Disord.* 2020;21(1):397.
- <sup>8</sup> National Joint Registry for England, Wales, Northern Ireland, the Isle of Man and Guernsey (NJR). Summary Report SP Humeral Affinis Inverse. 17.05.2023. Data valid to 12 May 2025.
- <sup>9</sup> NJR. National Joint Registry for England, Wales, Northern Ireland, the Isle of Man and Guernsey (NJR). Summary Report SP Humeral Affinis Inverse Cementless. 17.05.2023. Data valid to 12 May 2025.
- <sup>10</sup> Hamada K, Fukuda H, Mikasa M, Kobayashi Y. Roentgenographic findings in massive rotator cuff tears. A long-term observation. *Clin Orthop Relat Res.* 1990(254):92-6.
- <sup>11</sup> Smith PN, Gill DR, McAuliffe MJ, McDougall C, Stoney JD, Vertullo CJ, Wall CJ, Corfield S, Page R, Cuthbert AR, Du P, Harries D, Holder C, Lorimer MF, Cashman K, Lewis PL. Hip, Knee and Shoulder Arthroplasty: 2023 Annual Report, Australian Orthopaedic Association National Joint Replacement Registry, AOA: Adelaide, South Australia. 2023. <https://doi.org/10.25310/YWQZ9375>. Tables ST6 and ST56
- <sup>12</sup> Joudet T and Egger M. L'Inversion du couple frottement. In: L. Favard and P. Mansat, editors. *Prothèse totale d'épaule inversée.* Elsevier (SOFECOT): Elsevier; 2020.
- <sup>13</sup> <https://www.odep.org.uk/odep-products>, last access 16.08.2023.

AOANJRR-Tabellen <sup>11</sup>

**Table ST6 Cumulative Percent Revision of Primary Total Shoulder Replacement by Class (All Diagnoses)**

Shoulder Class	N Revised	N Total	1 Yr	3 Yrs	5 Yrs	7 Yrs	10 Yrs	14 Yrs
Total Stemmed Reverse	1728	47251	2.4 (2.2, 2.5)	3.5 (3.4, 3.7)	4.1 (3.9, 4.4)	4.7 (4.4, 4.9)	5.7 (5.3, 6.1)	7.3 (6.4, 8.3)

Note: Restricted to modern prostheses

**Table ST56 Cumulative Percent Revision of Primary Total Stemmed Reverse Shoulder Replacement by Prosthesis Combination**

Humeral Stem	Glenoid Component	N Revised	N Total	1 Yr	3 Yrs	5 Yrs	7 Yrs	10 Yrs	14 Yrs
Affinis	Affinis	62	1868	1.8 (1.3, 2.5)	3.5 (2.7, 4.6)	4.0 (3.1, 5.2)	5.0 (3.2, 7.6)		

Note: Restricted to modern prostheses

Only prostheses with >50 procedures have been listed

# *Preservation in motion*