


Manufacturer

Accuratus AG

Distributed by

MATHYS 
European Orthopaedics



For healthcare professional use only. The illustrated image does not represent a connection between the use of the medical device described, nor its performance.

Preservation in motion




Basic instruments hip

Basic instruments hip

Mathys offers various basic instruments to complement our hip portfolio. The bone retractors are available with **standard or long** soft tissue protection areas, and the **option of a handle** extension. The saddle-shaped geometry reduces trauma to the soft tissues. Due to the Subtilis 2.0 surface, the instruments are more corrosion-resistant, easier to clean and thus have a longer service life. The **second angle in the handle** also provide a better view of the operative field. The **subtle design** improves the approach compared to existing retractors, especially in narrow spaces. The **interconnecting radii** enable the best possible distribution of pressure on muscles and soft tissues.

In addition to retractors, we offer other basic instruments such as mallets, forceps, chisels, curettes, rasps and bone spreaders, which are helpful to the surgeon as basic instruments during treatment. Many of these instruments are modern optimisations of the designs made for Professor M. E. Müller in the 1950s and 1960s, developed in close cooperation with our clinical partners.



Bone retractors

- Pressure distribution on muscles and soft tissues optimised via saddle-shaped design with the option of a long soft tissue protection area and a handle extension
- Good view of the operative field due to glare-free Subtilis 2.0 surfaces and second angle in the handle
- Optimum access even in narrow spaces due to subtle design
- Pressure distribution on muscles and soft tissues optimised compared to conventional designs via two interconnecting radii

Advantages



Chisels

- High cutting performance and durability due to hardened, razor-sharp steel
- Good view of the operative field thanks to glare-free Subtilis 2.0 surfaces
- Durable, safe, hygienic and moisture-repellent plastic handles



Curettes and rasps


- Ergonomic modification offers better patient protection
- Good view of the operative field thanks to glare-free Subtilis 2.0 surfaces
- Durable, safe, hygienic and moisture-repellent plastic handles



Forceps and bone spreaders

- Lightweight yet durable design allows weight reduction in surgical trays
- Easy handling with intuitive design

Instruments



Lever

Item no.	Description
24.51.100	Extender Duplus 24 cm

Item no.	Description
24.51.37L	Lever with tissue protection 12x75/26 cm

Item no.	Description
24.51.38	Bone Lever Hohmann 18x40 mm/24 cm

Item no.	Description
24.51.38L	Lever with tissue protection 18x75/27 cm

Item no.	Description
24.51.40	Hip Lever 24x45 mm blunt/27 cm

Item no.	Description
24.51.40L	Hip Lever 24x80 mm blunt/32 cm

Item no.	Description
24.51.41	Bone Lever 24 mm blunt double-angled 25 cm

Item no.	Description
24.51.46	Pubis-Retractor, short beak 22/26 cm



Item no.	Description
24.51.48	Lever double angled 15x45 mm/28 cm



Item no.	Description
24.51.48L	Lever double angled 15x100 mm/31 cm



Item no.	Description
24.51.50	Bone Lever Amin narrow 22/29 cm



Item no.	Description
24.51.51	Bone Lever Amin large 32/30 cm



Item no.	Description
24.51.52	Bone Lever deep 22/28.5 cm



Item no.	Description
24.51.53	Bone lever 32/30 cm special Model



Item no.	Description
24.51.57	Bone Lever double angled 20mm/31 cm

Item no.	Description
24.51.58	Bone lever double angled 26 mm/31 cm

Item no.	Description
24.51.63L	Hohmann 41 x65 mm/32 cm

Item no.	Description
24.51.67	Soft tissue retractor 90° curved 40/30 cm

Item no.	Description
24.51.70	Troch. Lever 65°, V-spikes 30 mm/32 cm

Item no.	Description
24.51.71	Troch. Lever 75°, V-spikes 30 mm/30 cm

Item no.	Description
24.51.72	Troch. Lever 40°, V-spikes 30 mm/33 cm



Item no.	Description
24.51.73	Pelvis Lever 22 mm/24 cm



Item no.	Description
24.51.80	Soft Tissue Lever Nötzli left



Item no.	Description
24.51.81	Soft Tissue Lever Nötzli right




Item no.	Description
24.51.90	Pelvis Lever w/ plastic handle 26 mm/38 cm



Item no.	Description
24.51.91	Pelvis Lever 26 mm/38 cm w/ weight slot



Item no.	Description
36.60.03	Luxation lever 22 mm/31 cm



Item no.	Description
99.24.31	Hohmann 18x37 mm/24 cm/spec. curved



Bone preparation

Item no.	Description
16.20.51	Rongeur Stille-Luer curved 10 mm/22 cm

Item no.	Description
29.08.14-S	Rasp straight 14 mm round/20.5 cm

Item no.	Description
29.09.14-S	Rasp curved 14 mm round/20 cm

Item no.	Description
29.11.20-S	Rasp straight 20 mm round/20 cm

Item no.	Description
31.60.01	Bone Curette 9x12 mm/29 cm
31.60.02	Bone curette 15x20 mm/29 cm
31.60.11	Bone Curette toothed 9x12 mm/29 cm
31.60.12	Bone Curette toothed 15x20 mm/29 cm

Item no.	Description
34.00.35	Mallet with synthetic handle 300 gr.

Item no.	Description
34.00.55	Mallet with synthetic handle 500 gr.

Item no.	Description
34.00.75	Mallet with synthetic handle 750 gr.



Item no.	Description
35.12.08	Chisel Lexer w/ impact plate 8 mm/24.5 cm
35.12.15	Chisel Lexer straight 15 mm/23 cm

Item no.	Description
35.13.15	Chisel straight 16 mm/32 cm

Item no.	Description
35.13.16	Chisel curved 16 mm/32 cm

Item no.	Description
35.15.16	Gooseneck Chisel 16 mm/34.5 cm

Item no.	Description
35.15.28	Gooseneck Chisel 28 mm/32 cm

Item no.	Description
35.54.57	Osteotome Lambotte straight 6 mm/25 cm
35.54.66	Osteotome Lambotte straight 15 mm/25 cm
35.54.71	Osteotome Lambotte straight 20 mm/25 cm

Item no.	Description
36.50.10	Tamper metallic 17x7 mm concave/31 cm

Item no.	Description
36.50.11	Tamper metallic 17x7 mm straight/31 cm

Item no.	Description
36.60.02	Femoral Head Extractor, graduated/25 cm



Pliers & spreaders

Item no.	Description
24.01.58	Retractor Volkmann 1 tooth sharp 20 cm

Item no.	Description
24.01.61	Femurtoller angled 22 cm

Item no.	Description
24.01.73	Retractor Volkmann 4 teeth 9x19 mm/21.5 cm

Item no.	Description
24.02.68	Retractor Langenbeck 40x11 mm/21 cm

Item no.	Description
24.04.51-S	Retractor Kocher 40x25 mm/22.5 cm

Item no.	Description
24.04.52-S	Retractor Kocher 70x25 mm/24 cm

Item no.	Description
24.04.53-S	Retractor Kocher 50x25 mm/24 cm



Item no.	Description
27.10.16	Self-retaining Retractor 3x4 teeth, 16 cm

Item no.	Description
27.20.51	Self-retaining Retractor 4x4 teeth, 31 cm

Dissection

Item no.	Description
02.00.03	Scalpel Handle Nr. 3

Item no.	Description
02.00.04	Scalpel Handle Nr. 4

Item no.	Description
02.01.03	Scalpel Handle Nr. 3L

Item no.	Description
02.01.04	Scalpel Handle Nr. 4L



Item no.	Description
04.15.14	Op Scissors Mayo straight blunt 14 cm

Item no.	Description
04.16.14	Op Scissors Mayo curved blunt 14 cm



Item no.	Description
04.30.14	Op Scissors Metzenbaum curved blunt 14 cm

Item no.	Description
04.30.20	Op Scissors Metzenbaum curved blunt 20 cm

Item no.	Description
04.90.95	Scissors curved, blunt, 23 cm / mod.Nötzli

Item no.	Description
06.00.14	Dis Forceps standard straight 14 cm

Item no.	Description
06.00.20	Dis Forceps standard straight 20 cm

Item no.	Description
07.60.15	Tissue Forceps Gillies straight 15 cm

Item no.	Description
20.00.16	Needle Holder w/ Inserts Mayo-Hegar 16 cm

Item no.	Description
20.00.20	Needle Holder w/ Inserts Mayo-Hegar 20 cm

Manufacturer

Accuratus AG

Distributed by

MATHYS 
European Orthopaedics

Australia	Mathys Orthopaedics Pty Ltd Artarmon, NSW 2064 Tel: +61 2 9417 9200 info.au@mathysmedical.com	Italy	Mathys Ortopedia S.r.l. 20141 Milan Tel: +39 02 4959 8085 info.it@mathysmedical.com
Austria	Mathys Orthopädie GmbH 2351 Wiener Neudorf Tel: +43 2236 860 999 info.at@mathysmedical.com	Japan	Mathys KK Tokyo 108-0075 Tel: +81 3 3474 6900 info.jp@mathysmedical.com
Belgium	Mathys Orthopaedics Belux N.V.-S.A. 3001 Leuven Tel: +32 16 38 81 20 info.be@mathysmedical.com	New Zealand	Mathys Ltd. Auckland Tel: +64 9 478 39 00 info.nz@mathysmedical.com
France	Mathys Orthopédie S.A.S 63360 Gerzat Tel: +33 4 73 23 95 95 info.fr@mathysmedical.com	Netherlands	Mathys Orthopaedics B.V. 3001 Leuven Tel: +31 88 1300 500 info.nl@mathysmedical.com
Germany	Mathys Orthopädie GmbH «Centre of Excellence Sales» Bochum 44809 Bochum Tel: +49 234 588 59 0 sales.de@mathysmedical.com «Centre of Excellence Ceramics» Mörsdorf 07646 Mörsdorf/Thür. Tel: +49 364 284 94 0 info.de@mathysmedical.com «Centre of Excellence Production» Hermsdorf 07629 Hermsdorf Tel: +49 364 284 94 110 info.de@mathysmedical.com	P. R. China	Mathys (Shanghai) Medical Device Trading Co., Ltd Shanghai, 200041 Tel: +86 21 6170 2655 info.cn@mathysmedical.com
		Switzerland	Mathys (Schweiz) GmbH 2544 Bettlach Tel: +41 32 644 1 458 info@mathysmedical.com
		United Kingdom	Mathys Orthopaedics Ltd Alton, Hampshire GU34 2QL Tel: +44 8450 580 938 info.uk@mathysmedical.com

Local Marketing Partners in over 30 countries worldwide ...

Manufactured by Accuratus AG • Holligenstrasse 94 • 3008 Bern • Switzerland • www.accuratus.ch

Distributed by Mathys Ltd Bettlach • Robert Mathys Strasse 5 • P.O. Box • 2544 Bettlach • Switzerland • www.mathysmedical.com