

Manufacturer

Effectum Medical AG


Distributed by

MATHYS 
European Orthopaedics



For healthcare professional use only. The illustrated image does not represent a connection between the use of the medical device described, nor its performance.

Preservation in motion




RAP-hip

RAP-hip

RAP hip instruments enable revision of a wide range of hip stems using **universal and special adapters** and can be used via various approaches. Surgery time can be reduced thanks to its **easy application and simple assembly / disassembly**. The surgeon can generate a **high axial force** with the integrated guided slide weight, resulting in an **excellent separation effect between bone and implant**.

Advantages



- Revision of a wide range of hip stems with use of universal and special adapters
- Can be inserted via dorsal, lateral or anterior approaches
- Reduced surgery time due to easy application and simple assembly/disassembly



- High axial force due to integrated guided slide weight
- Avoids power loss and vibrations at the wrist joint thanks to the patented handle damping system

Instruments



Item no.	Description
300-00116	Extension rod
300-00118	M8 screws
300-00122	Large impact piece 860 g
300-00180	Handle for guide rod
300-00185	Nut
350-00145	Impactor
350-00183	Guide rod
350-00193	Compression spring with washer
350-00150	Tool for basic set
500-00191	11 mm AF open-end/ring wrench



Item no.	Description
350-00162	Femoral head extractor fork
300-00164	Femoral head extractor mandrel
300-00236	Femoral head extractor fixation



Item no.	Description
350-00222	M6 special adapter
350-00226	M8 special adapter
350-00375	M7 special adapter
350-00377	1/4"-20UNC special adapter
330-00340	Insertion aid



Item no.	Description
350-00148	Universal adapter module 12/14
350-00152	Universal adapter module 14/16
350-00156	Universal adapter module 10/12
350-00344	Universal adapter module 8/10
300-00137	Pressure screw
300-00138	Universal adapter mandrel
300-00190	Universal adapter handle

Item no.	Description
300-00314	Short hook
300-00315	PROTEK connector

Manufacturer

Effectum Medical AG

Distributed by



Australia Mathys Orthopaedics Pty Ltd
Artarmon, NSW 2064
Tel: +61 2 9417 9200
info.au@mathysmedical.com

Austria Mathys Orthopädie GmbH
2351 Wiener Neudorf
Tel: +43 2236 860 999
info.at@mathysmedical.com

Belgium Mathys Orthopaedics Belux N.V.-S.A.
3001 Leuven
Tel: +32 16 38 81 20
info.be@mathysmedical.com

France Mathys Orthopédie S.A.S
63360 Gerzat
Tel: +33 4 73 23 95 95
info.fr@mathysmedical.com

Germany Mathys Orthopädie GmbH
«Centre of Excellence Sales» Bochum
44809 Bochum
Tel: +49 234 588 59 0
sales.de@mathysmedical.com

«Centre of Excellence Ceramics» Mörsdorf
07646 Mörsdorf/Thür.
Tel: +49 364 284 94 0
info.de@mathysmedical.com

«Centre of Excellence Production» Hermsdorf
07629 Hermsdorf
Tel: +49 364 284 94 110
info.de@mathysmedical.com

Italy Mathys Ortopedia S.r.l.
20141 Milan
Tel: +39 02 4959 8085
info.it@mathysmedical.com

Japan Mathys KK
Tokyo 108-0075
Tel: +81 3 3474 6900
info.jp@mathysmedical.com

New Zealand Mathys Ltd.
Auckland
Tel: +64 9 478 39 00
info.nz@mathysmedical.com

Netherlands Mathys Orthopaedics B.V.
3001 Leuven
Tel: +31 88 1300 500
info.nl@mathysmedical.com

P. R. China Mathys (Shanghai) Medical Device Trading Co., Ltd
Shanghai, 200041
Tel: +86 21 6170 2655
info.cn@mathysmedical.com

Switzerland Mathys (Schweiz) GmbH
2544 Bettlach
Tel: +41 32 644 1 458
info@mathysmedical.com

United Kingdom Mathys Orthopaedics Ltd
Alton, Hampshire GU34 2QL
Tel: +44 8450 580 938
info.uk@mathysmedical.com

Local Marketing Partners in over 30 countries worldwide ...

Manufactured by Effectum Medical AG • Kirchgasse 11 • 4600 Olten • Switzerland • www.effectummedical.com

Distributed by Mathys Ltd Bettlach • Robert Mathys Strasse 5 • P.O. Box • 2544 Bettlach • Switzerland • www.mathysmedical.com