

Manufacturer

Dedienne Santé

Distributed by

MATHYS 
European Orthopaedics

New Dynamic Trial Inlays for DS Evolution

The DS Evolution system has new Dynamic Trial Inlays.



Trial liners are available to test the joint's stability and range of motion.



Dynamic trials can be carried out after the metal cup has been impacted.



Item no.	Description	Size
3700502204674	Trial Liner 46/28	46
3700502204681	Trial Liner 48/28	48
3700502204698	Trial Liner 50/28	50
3700502204704	Trial Liner 52/28	52
3700502204711	Trial Liner 54/28	54
3700502204728	Trial Liner 56/28	56
3700502204735	Trial Liner 58/28	58
3700502204742	Trial Liner 60/28	60
3700502204759	Trial Liner 62/28	62
3700502204766	Trial Liner 64/28	64