

For healthcare professional use only. The illustrated image does not represent a connection between the use of the medical device described, nor its performance.



balanSys REV


balanSys REV

Knee revision system

Proven articulation

balanSys REV Femur


- Identical articulation as balanSys BICONDYLAR
- **Compatible** to balanSys PS Inlays for complex primary TKA*
- **Compatible** to all balanSys Patellas
- 5 sizes A/B/C/D/E
- Left and right



Stability

balanSys REV Inlays

- Up to $\pm 4^\circ$ rotational freedom
- Post to shaft connection **reinforced by metal**
- **Jumping height** 21 mm
- Composition up to 33 mm (incl. augment)
- 6 heights 10.5/13/15.5/18/20.5/23



Variability

balanSys REV Stem

- **Identical** stems for tibia and femur
- With morse taper connection
- **Straight** and 4 mm **offset**
- 12 diameters: 10 to 24mm
- 3 lengths: 80, 140 and 200 mm

* For combination possibilities to balanSys BICONDYLAR Implants refer to surgical technique balanSys REV, Item no. 336.030.052 or 336.030.049



Bone defect management
balanSys REV Femur Augmentation

- Distal and dorsal, 5 and 10 mm

Proven geometry
balanSys REV Tibia Plateau

- **Compatible** to balanSys PS, CR, UC Inlays for complex primary TKA*
- 5 sizes 64/70/75/80/85

Bone defect management
balanSys REV Tibia Augmentation

- Half size blocks with **8° inclination**
- 5 and 10 mm

Precision
balanSys REV Instruments

- Reproducible orientation of stems
- All cuts guided

Australia	Mathys Orthopaedics Pty Ltd Lane Cove West, NSW 2066 Tel: +61 2 9417 9200 info.au@mathysmedical.com	Italy	Mathys Ortopedia S.r.l. 20141 Milan Tel: +39 02 5354 2305 info.it@mathysmedical.com
Austria	Mathys Orthopädie GmbH 2351 Wiener Neudorf Tel: +43 2236 860 999 info.at@mathysmedical.com	Japan	Mathys KK Tokyo 108-0075 Tel: +81 3 3474 6900 info.jp@mathysmedical.com
Belgium	Mathys Orthopaedics Belux N.V.-S.A. 3001 Leuven Tel: +32 16 38 81 20 info.be@mathysmedical.com	New Zealand	Mathys Ltd. Auckland Tel: +64 9 478 39 00 info.nz@mathysmedical.com
France	Mathys Orthopédie S.A.S 63360 Gerzat Tel: +33 4 73 23 95 95 info.fr@mathysmedical.com	Netherlands	Mathys Orthopaedics B.V. 3001 Leuven Tel: +31 88 1300 500 info.nl@mathysmedical.com
Germany	Mathys Orthopädie GmbH «Centre of Excellence Sales» Bochum 44809 Bochum Tel: +49 234 588 59 0 sales.de@mathysmedical.com «Centre of Excellence Ceramics» Mörsdorf 07646 Mörsdorf/Thür. Tel: +49 364 284 94 0 info.de@mathysmedical.com «Centre of Excellence Production» Hermsdorf 07629 Hermsdorf Tel: +49 364 284 94 110 info.de@mathysmedical.com	P. R. China	Mathys (Shanghai) Medical Device Trading Co., Ltd Shanghai, 200041 Tel: +86 21 6170 2655 info.cn@mathysmedical.com
		Switzerland	Mathys (Schweiz) GmbH 2544 Bettlach Tel: +41 32 644 1 458 info@mathysmedical.com
		United Kingdom	Mathys Orthopaedics Ltd Alton, Hampshire GU34 2QL Tel: +44 8450 580 938 info.uk@mathysmedical.com

Local Marketing Partners in over 30 countries worldwide ...

